

High-quality Nucleic Acids from Pathogen-infected Stool Samples with GENTi™ Advanced Viral DNA/RNA Extraction Kit

Experimental Conditions

Materials Required

- ♦ GENTi™³² Advanced Automatic Extraction System (GTI032A)
- ♦ GENTi™ Advanced Viral DNA/RNA Extraction Kit (902-096A)
- ♦ Stainless bead, 6 mm (SUS 304, Supplier : S)
- ♦ 1X PBS (Phosphate-buffered saline)
- ♦ 5 ml conical tube
- ♦ Pipette & sterile pipette tips
- ♦ Suitable protector (e.g., lab coat, disposable gloves, goggles, etc.)

Sample Information

- ♦ Sample type : Stool infected with pathogen
 - Viral RNA : Infectious bronchitis virus (IBV)
 - Bacterial DNA : *Mycoplasma gallisepticum* (MG)
- ♦ Extraction conditions
 - Sample amount : 200 µl
 - Elution volume : 100 µl
 - Extraction protocol : Viral_Normal (operation time : 29' 35")

Sample Preparation

1. Prepare each 200 mg stool infected with IBV or MG and transfer the samples to each 5 ml conical tube with 1X PBS and 6 mm stainless bead.
2. Vortex for 1 min or until the stool sample is thoroughly homogenized and incubate for 1 min at room temperature.
3. The next step is according to **GENTi™ Advanced Viral DNA/RNA Extraction Kit protocol**.

Protocol

GENTi™ Advanced Viral DNA/RNA Extraction Kit protocol

* For more details and methods, please refer to [the handbook of GENTi™ Advanced Viral DNA/RNA Extraction Kit](#).

1. Add 7 µl of dissolved Carrier RNA (1 µg/µl) to 1st/7th well of plate type cartridge.
2. Add 200 µl of samples to 1st/7th well.
3. Load the plate type cartridge on the tray of GENTi™³² Advanced Automatic Extraction System.
4. Insert the magnetic rod cover to the end to strip bracket.
5. Select the correct extraction protocol and operate the extraction system.

Result




Figure 1. Result of C_q value of DNA template extracted from infected stool.

The DNA templates were extracted from the stool infected with *Mycoplasma gallisepticum* (MG) using GENTi™ Advanced Viral DNA/RNA Extraction Kit on GENTi™³² Advanced Automatic Extraction System. Eluted DNA templates were analyzed with a TaqMan-based real-time PCR assay using CFX-96. Each eluate was analyzed in duplicate.

- PCR instrument : CFX-96 (1855201)
- qPCR kit : Probe qPCR Mix (RR391A)
- Target gene : None specific




Figure 2. Result of C_q value of RNA template extracted from infected stool.

The RNA templates were extracted from the stool infected with Infectious bronchitis virus (IBV) using GENTi™ Advanced Viral DNA/RNA Extraction Kit on GENTi™³² Advanced Automatic Extraction System. Eluted RNA templates were synthesized to cDNA with reverse transcription; and then analyzed with TaqMan-based one-step RT-qPCR assay using CFX-96. Each eluate was analyzed in duplicate.

- PCR instrument : CFX-96 (1855201)
- RT-qPCR kit : 2X 1 Step RT-qPCR Master Mix (for probe) (QRT1-XV-100R)
- Target gene : None specific