



Customer & Technical Support

Do not hesitate to ask us any question.

We thank you for any comment or advice.

Contact us at

www.geneall.com

Tel: 82-2-407-0096

Fax: 82-2-407-0779

E-mail(Order/Sales) : sales@geneall.com

E-mail(Tech. Info.): tech@geneall.com

Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

Trademarks

Hybrid- Q^{TM} , ExprepTM, ExfectionTM, ExpinTM, ExgeneTM, GenExTM, DirExTM, Hybrid-RTM, RiboExTM, RiboSaverTM, RiboSaverTM, EzClearTM, EzPureTM, EzPureTM, EzPassTM, AmpONETM, AmpMasterTM, HyperScriptTM, ProtinExTM, PAGESTATM, STEADiTM are trademarks of GeneAll Biotechnology co., Itd.

© 2016 GeneAll Biotechnology, all right reserved.

This protocol handbook is included in:

GeneAll® Exgene[™] FFPE Tissue DNA (138-150, 135-152)

Visit www.geneall.com or www.geneall.co.kr for FAQ, QnA and more information.

INDEX

Kit Contents	04
Quality Control	05
Storage Condition	
User Precautions	
Preparation of Proteinase K solution	06
Materials Not Provided	
Product Description	07
Protocol	80
Total DNA from FFPE Tissue	
Troubleshooting Guide	- 11
APPENDIX	13
Protocol for total DNA from FFPE Tissue using Xylene	
Ordering Information	17

Kit Contents

Exgene™ FFPE Tissue DNA

Components	Qua	ıntity	Storage
Components	138-150	138-152	Storage
Buffer DP	55 ml	140 ml x 2	
Buffer FPL	15 ml	60 ml	
Buffer FPB	15 ml	60 ml	
Buffer BW (concentrate) *	16 ml	90 ml	
Buffer TW (concentrate) *	10 ml	40 ml	
Buffer AE **	15 ml	60 ml	Room temperature
RNase A solution (100 mg/ml) ***	300 ul	1.2 ml	$(15 \sim 25^{\circ}C)$
Proteinase K (mg) ****	24 mg	120 mg	
PK Storage buffer	1.5 ml	7 ml	
GeneAll® Column type G (with collection tube)	50	250	
Collection tube	150	750	
Protocol Handbook	1	1	

^{*} Before using for the first time, add absolute ethanol (ACS grade or better) into buffer BW and TW as indicated on the bottle

Product Specifications

Exgene [™] FFPE Tissue DNA						
Туре	Spin					
Maximum amount of starting samples	Up to 8 sections of 10 um in thickness					
Maximum loading volume	~ 750 ul					
Minimum elution volume	~ 30 ul					
Maximum binding capacity	\sim 60 ug					

^{*} Contains sodium azide as a preservative

^{** 10} mM Tris-HCl, pH 9.0, 0.5 mM EDTA

^{***} For the long-term storage of RNase A, store below $4^{\circ}C$

^{****} After reconstitution of Proteinase K with PK Storage buffer, store at 4° C for its stability But for long-term storage, store at -20° C

Quality Control

All components in GeneAll[®] Exgene™ FFPE Tissue DNA kit are manufactured in strictly clean condition, and its degree of cleanness is monitored periodically. For consistency of product, the quality certification process is carried out thoroughly and only the qualified is delivered.

Storage Conditions

All components of GeneAll® Exgene TFPE Tissue DNA kit should be stored at room temperature (15 \sim 25°C). RNase A solution is delivered under ambient conditions and can be stored at room temperature for 6 months without significant decrease in activity. But for the long-term storage, storage at -20 \sim 8°C is recommended. After reconstitution of Proteinase K with storage buffer, it should be stored below 4°C for conservation of activity. It can be stored at 4°C for 1 year without significant decrease in activity. But for the long-term storage, storage at below -20°C is recommended.

During shipment or storage below cold ambient condition, a precipitate can be formed in buffer FPL or FPB. In such a case, heat the bottle at 56°C to dissolve thoroughly. Using precipitated buffers will lead to poor DNA recovery.

GeneAll® Exgene™ FFPE Tissue DNA kit is guaranteed for 18 months from the production date.

User Precautions

Before use, read these instructions carefully.

The buffers included in GeneAll® Exgene™ FFPE Tissue DNA kit contain irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling. Always wear gloves and eye protector, and follow standard safety precautions. In case of contact, wash immediately with plenty of water and seek medical advice.

In detail, Buffer DP contains irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling.

Buffer FPB and BW contain chaotropes. It can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste. Please consult the Material Safety Data Sheet (MSDS) for information regarding hazards and safe handling practices.

Preparation of Proteinase K solution (20 mg/ml) at first use

Proteinase K is provided in a freeze-dried format. Thus, it should be reconstituted thoroughly with PK Storage buffer before experiment. PK storage buffer contains calcium chloride and glycerol which do not affect the enzyme activity, but contribute to its stability.

To obtain a solution of 20 mg/ml of Proteinase K, add 1.2 ml (Cat. No.138-150) or 6 ml (Cat. No.138-152) of the PK Storage buffer to the glass vial containing lyophilized Proteinase K, and mix carefully and gently to avoid foaming.

Reconstituted enzyme should be stored at 4° C for its stability. But for long-term storage, storage at -20° C is recommended.

Materials Not Provided

Reagent

· Absolute ethanol, ACS grade or better

Disposable material

- Pipette tips
- 1.5 ml or 2.0 ml microcentrifuge tubes

Equipment

- Microtome
- Microcentrifuge
- Heating block or water bath
- Suitable protector (ex; lab coat, disposable gloves, goggles, etc)

Product Description

GeneAll[®] Exgene™ FFPE Tissue DNA kit provides a convenient and easy method for the isolation of total DNA from Formalin Fixed and Paraffin Embedded (FFPE) specimen by non-organic solvent. FFPE is one of the most commonly used methods of clinical tissue preservation; the clinical tissue is fixed by formalin and subsequently embedded in paraffin to keep its original form.

The FFPE tissue is useful in disease research such as microscopic observation and immunohistochemical analysis. And the extracted-nucleic acid from FFPE specimen can be used for molecular diagnosis of various diseases. However, during the fixative process, the nucleic acids in FFPE are damaged significantly by various degrees of crosslinking between DNA and protein, and the damage get worse during its long-term preservation. For such a reason, the DNA isolated from the preserved FFPE specimen generally has low qualities in its yield, purity, integrity and PCR-processivity. But despite these problems, the purified nucleic acids from FFPE specimen are widely used for the PCR targeted to relatively short DNA fragments.

To obtain DNA from FFPE tissue by GeneAll® Exgene™ FFPE Tissue DNA kit, FFPE specimen is deparaffinized in buffer DP which rapidly separate tissue from paraffin sections, and then the sample is lysed in the optimized buffer containing detergents and lytic enzymes. Under high salt condition, DNA in the lysate binds to silica membrane and impurities pass through membrane in to a collection tube. The membrane is washed with a series of alcohol-containing buffer to remove any traces of proteins, cellular debris and salts. Finally pure DNA is released into a clean collection tube with deionized water or low ionic strength buffer.

This fast and convenient procedure of GeneAll[®] Exgene[™] FFPE Tissue DNA kit takes only 150 min for complete preparations of total DNA. Purified DNA can be used directly for PCR (\leq 500 bp), real-time PCR, and other downstream applications.

GeneAll® Protocol for FFPE Tissue

Before experiment

- Prepare the heat block or water bath to 56 °C and 90 °C
- Prepare absolute ethanol
- Prepare 1.5 ml and 2.0 ml microcentrifuge tube
- All centrifugation should be performed at room temperature
- Buffer FPL and FPB may precipitate at cool ambient temperature If so, dissolve it in 56 °C water bath
- 1. Cut up to 8 sections of 5 10 um in thickness from the FFPE sample block.

The amount of starting sample should not exceed 8 sections of 10 um in thickness. It is recommended to use $1 \sim 3$ sections of $5 \sim 10$ um in thickness.

2. Place the sections in 2.0 ml microcentrifuge tube.

Trim the marginal paraffin off from the tissue sections as much as possible.

- 3. Add 1000 ul of buffer DP and vortex to mix.
- 4. Incubate at 56°C for 3 min and spin down briefly to remove any drops from inside of the lid.

Make sure that paraffin thoroughly melts during the incubation step.

5. Carefully discard the buffer DP as much as possible by pipetting without any loss of tissue pieces.

Too much paraffin in the starting sample may cause re-solidification of melted paraffin. If so, repeat step $3\sim 5$ once more to remove the residual paraffin from the sample.

6. Add 180 ul of buffer FPL and mix thoroughly by vigorous vortexing.

7. Add 20 ul of Proteinase K solution (20 mg/ml, provided) and mix thoroughly by vortexing or pipetting. Incubate at 56°C for 1 hr.

It is essential to mix the components thoroughly for proper lysis.

Lysis time varies from 1 hr to overnight depending on the type of tissue and the starting amount. The lysate should become translucent without any particles after complete lysis. To enhance the lysis efficiency during this incubation step, vortex the tube every 15 min or use shaking incubator or agitator.

8. Incubate at 90°C for I hr and spin down briefly to remove any drops from inside of the lid.

This step allows DNA to be decrosslinked from DNA-Protein crosslinking. Cool down to room temperature before proceeding to next step.

- (Optional:) If RNA-free DNA is required, add 4 ul of RNase A solution (100 mg/ml, provided), vortex to mix thoroughly, and incubate for 2 min at room temperature.
- 10. Add 200 ul of buffer FPB to the tube and mix thoroughly by vortexing. Spin down briefly to remove any drops from inside of the lid.

Disregard a little residual buffer DP in upper phase, because it will not affect the next procedure. If the lysate volume is larger than 200 ul, the volume of buffer FPB should be adjusted proportionally.

II. Add 200 ul of absolute ethanol (not provided) to the sample, pulse-vortex to mix the sample thoroughly, and spin down briefly to remove any drops from inside of the lid.

If the lysate volume is larger than 200 ul, the volume of ethanol should be adjusted proportionally.

12. Transfer all of the mixture to the SV column carefully, centrifuge for 1 min at 6,000 x g above (>8,000 rpm) and replace the collection tube with new one (provided).

If the mixture has not passed completely through the membrane, centrifuge again at full speed ($> 13,000 \times g$) until all of the solution has passed through. Centrifugation at full speed will not adversely affect DNA recovery.

13. Add 600 ul of buffer BW, centrifuge for I min at 6,000 x g above (>8,000 rpm) and replace the collection tube with new one (provided).

If the SV column has a colored residue after centrifugation, repeat this step until no colored residue remains. See Trouble shooting guide for detail.

Centrifugation at full speed will not adversely affect DNA recovery.

14. Add 700 ul of buffer TW. Centrifuge for I min at 6,000 x g above (>8,000 rpm). Discard the pass-through and reinsert the SV column back into the collection tube.

Centrifugation at full speed will not adversely affect DNA recovery.

15. Centrifuge at full speed (>13,000 x g) for I min to remove residual wash buffer. Place the SV column into a fresh 1.5 ml microcentrifuge tube (not provided).

Care must be taken at this step for eliminating the carryover of buffer TW. If a carryover of buffer TW still occurs, centrifuge again for I min at full speed with the collection tube before transferring to a new 1.5 ml microcentrifuge tube.

Centrifugation must be performed at full speed (13,000 x g \sim 20,000 x g).

16. Add 50 ul of buffer AE or distilled deionized water. Incubate for 1 min at room temperature. Centrifuge at full speed (>13,000 x g) for 1 min.

Elution volume can be adjusted according to an experiment's purpose.

For long-term storage of purified DNA, eluting in buffer AE is recommended. But EDTA included in buffer AE can inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled deionized water (>pH7.0) or Tris-HCl (>pH8.5).

When using water for elution, check the pH of water before elution.

Troubleshooting Guide

Facts	Possible Causes	Suggestions
Low or no recovery	Insufficient lysis	Incomplete lysis is due to too much starting material, imperfect mixing with buffer FPL, insufficient time for complete lysis or not enough disruption of sample.
	Weakened activity of Proteinase K caused by inappropriate storage condition or out-of date	Proteinase K must be stored below 4°C for maintenance of proper activity. Lysis cannot be done properly with degenerated Proteinase K. Replace with new one.
	Starting sample is too old or inappropriately fixed	Too old or inappropriately fixed sample often cause low recovery. If possible, use the FFPE sample fixed within I hour of surgical resection and use fresh FFPE sample.
	Loss of starting sample during deparaffinization	Intensive care must be taken when removing the buffer DP with paraffin in order not to lose the starting sample.
SV column has colored residue associated with it after wash, resulting in colored residue	Insufficient lysis	Insufficient lysis may cause not complete removal of colored residue on SV column membrane. Repeat the procedure after refer to trouble shooting of 'insufficient lysis' of "Low or no recovery".
	Incomplete removal of hemoglobin	In case of DNA preparation from certain animal blood, it is hard to remove hemoglobin from their blood. Carry out additional wash step with buffer BW.
Column clogging	Sample not homogenized completely	Confirm the complete lysis of sample in buffer FPL. And make sure to incubate the lysate with occasional vortexing.
	Too much starting sample	Reduce the amount of starting sample.

Facts	Possible Causes	Suggestions
Low A _{260/280}	Sample not homogenized completely	Confirm the complete lysis of sample in buffer FPL. And make sure to incubate the lysate with occasional vortexing.
Too much degraded DNA	Starting sample is too old or inappropriately fixed	Using too old sample or inappropriately fixed sample often results in extreme degradation of purified DNA. If possible, use the fresh sample which fixed within an hour after a surgical resection.
PCR reaction is not performed well with purified DNA	DNA fragmentation	Using too old sample or inappropriately fixed sample often results in extreme degradation of purified DNA. It is strongly recommended to design PCR primers to target the fragment as short as possible (< 500 nucleotides). And if possible, use the fresh FFPE sample fixed within 1 hour after a surgical resection.
	Incompletely decrosslinked	FFPE samples are crosslinked between formalin-DNA-protein. The purified crosslinked-DNA can attenuate the processivity of PCR reaction. It is strongly recommended to design PCR primers to target the fragment as short as possible (< 500 nucleotides).

APPENDIX

Protocol for total DNA from FFPE Tissue using Xylene

Before experiment

- Prepare the heat block or water bath to 56 °C and 90 °C
- Prepare xylene
- Prepare absolute ethanol
- Prepare 1.5 ml and 2.0 ml microcentrifuge tubes
- All centrifugation should be performed at room temperature
- Buffer FPL and FPB may precipitate at cool ambient temperature
 If so, dissolve it in 56°C water bath
- 1. Cut up to 8 sections of 5 10 um in thickness from the FFPE sample block.

The amount of starting sample should not exceed 8 sections of 10 um in thickness. It is recommended to use $1\sim3$ sections of $5\sim10$ um in thickness.

2. Place the sections in 2.0 ml microcentrifuge tube.

Trim the marginal paraffin off from the tissue sections as much as possible.

- 3. Add 1200 ul xylene and vortex to mix.
- 4. Incubate at 56°C for 10 min and centrifuge at full speed (>13,000 x g) for 5 min. Carefully discard supernatant by pipetting.

Make sure that paraffin thoroughly melts during the incubation step. Be careful not to lose any pellet.

- 5. Add 1200 ul of absolute ethanol to the pellet to discard residual xylene and mix by vortexing.
- Centrifuge at full speed for 5 min and carefully discard the ethanol by pipetting.

Do not remove any tissues.

- 7. Repeat the steps $5 \sim 6$ once.
- 8. Open the tube and incubate at room temperature or 37° C for $10 \sim 15$ min.

Residual ethanol may interfere with purification of nucleic acid. Make sure that ethanol in tube has to be evaporated thoroughly.

- 9. Add 180 ul of buffer FPL and mix thoroughly by vigorous vortexing.
- 10. Add 20 ul of Proteinase K solution (20 mg/ml, provided) and mix thoroughly by vortexing or pipetting. Incubate at 56°C for 1 hr.

It is essential to mix the components thoroughly for proper lysis.

Lysis time varies from 1 hr to overnight depending on the type of tissue and the starting amount. The lysate should become translucent without any particles after complete lysis. To enhance the lysis efficiency during this incubation step, vortex the tube every 15 min or use shaking incubator or agitator.

11. Incubate at 90°C for 1 hr and spin down briefly to remove any drops from inside of the lid.

This step allows DNA to be decrosslinked from DNA-Protein crosslinking. Cool down to room temperature before proceeding to next step.

- 12. (Optional:) If RNA-free DNA is required, add 4 ul of RNase A solution (100 mg/ml, provided), vortex to mix thoroughly, and incubate for 2 min at room temperature.
- 13. Add 200 ul of buffer FPB to the tube and mix thoroughly by vortexing. Spin down briefly to remove any drops from inside of the lid.

Disregard a little residual buffer DP in upper phase, because it will not affect the next procedures. If the lysate volume is larger than 200 ul, the volume of buffer FPB can be adjusted proportionally.

14. Add 200 ul of absolute ethanol (not provided) to the sample, pulse-vortex to mix the sample thoroughly, and spin down briefly to remove any drops from inside of the lid.

If the lysate volume is larger than 200 ul, the volume of ethanol should be adjusted proportionally.

15. Transfer all of the mixture to the SV column carefully, centrifuge for I min at 6,000 x g above (>8,000 rpm) and replace the collection tube with new one (provided).

If the mixture has not passed completely through the membrane, centrifuge again at full speed ($>13,000 \times g$) until all of the solution has passed through. Centrifugation at full speed will not adversely affect DNA recovery.

16. Add 600 ul of buffer BW, centrifuge for 1 min at 6,000 x g above (>8,000 rpm) and replace the collection tube with new one (provided).

If the SV column has a colored residue after centrifugation, repeat this step until no colored residue remains. See Trouble shooting guide for detail.

Centrifugation at full speed will not adversely affect DNA recovery.

17. Add 700 ul of buffer TW. Centrifuge for 1 min at 6,000 x g above (>8,000 rpm). Discard the pass-through and reinsert the SV column back into the collection tube.

Centrifugation at full speed will not adversely affect DNA recovery.

18. Centrifuge at full speed (>13,000 x g) for 1 min to remove residual wash buffer. Place the SV column into a fresh 1.5 ml microcentrifuge tube (not provided).

Care must be taken at this step for eliminating the carryover of buffer TW. If a carryover of buffer TW still occurs, centrifuge again for I min at full speed with the collection tube before transferring to a new I.5 ml microcentrifuge tube.

Centrifugation must be performed at full speed (13,000 x g \sim 20,000 x g).

19. Add 50 ul of buffer AE or distilled deionized water. Incubate for I min at room temperature. Centrifuge at full speed (>13,000 x g) for 1 min.

Elution volume can be adjusted according to an experiment's purpose.

For long-term storage of purified DNA, eluting in buffer AE is recommended. But EDTA included in buffer AE can inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled deionized water (>pH7.0) or Tris-HCI (>pH8.5). When using water for elution, check the pH of water before elution.

Ordering Information

Products	Scale	Size	Cat. No.	Туре	Products	Scale	Size	Cat. No.	Туре
				71					71
GeneAll® Hybri o	1-Q™ fo	r rapid p	reparation of	blasmid DNA	GeneAll® Exgene	for is	olation o	f total DNA	
Plasmid Rapidprep	mini	50	100-150	spin		mini	100	105-101	spin /
		200	100-102		-		250	105-152	vacuum
C A II® F				2010	Blood SV	Midi	26	105-226	spin /
GeneAll® <i>Expre</i>	p for pi	reparatio 50		JNA	-		100	105-201	vacuum · ,
	mini	200	101-150	spin /		MAXI	26	105-316	spin / vacuum
	1111111	1,000	101-111	vacuum			100	106-101	spin /
Plasmid SV		26	101-226			mini	250	106-152	vacuum
	Midi	50	101-250	spin /	Cell SV -	N44371	10	106-310	spin /
		100	101-201	vacuum		MAXI	26	106-326	vacuum
						mini	100	108-101	spin /
GeneAll® <i>Exfect</i>		tuan-f-	tion and to the	anaid DAIA	_	1111111	250	108-152	vacuum
for prepa	aration of		tion-grade pla		Clinic SV	Midi	26	108-226	spin /
Diameted I E	mini	50 200	111-150	spin / vacuum			100	108-201	vacuum
Plasmid LE (Low Endotoxin)		26	111-102			MAXI	10	108-310	spin /
,	Midi	100	111-220	spin / vacuum			26	108-326	vacuum
Plasmid EF		20	121-220		Genomic DNA micro)	100	118-050	spin
(Endotoxin Free)	Midi	100	121-201	spin		mini	250	117-101	spin / vacuum
				-		26	117-132	spin /	
GeneAll® <i>Expin</i> ¹	or pur	ification	of fragment D	NA	Plant SV	Midi	100	117-220	vacuum
C-10/		50	102-150	spin /	-		10	117-310	spin /
Gel SV	mini	200	102-102	vacuum		MAXI	26	117-326	vacuum
PCR SV	mini	50	103-150	spin /	Soil DNA mini	mini	50	114-150	spin
	1111111	200	103-102	vacuum	Stool DNA mini	mini	50	115-150	spin
CleanUp SV	mini	50	113-150	spin /	Viral DNA / RNA	mini	50	128-150	spin
		200	113-102	vacuum	FFPE tissue DNA	mini	50	138-150	spin
Combo GP	mini	50	112-150	spin /			250	138-152	Sp
		200	112-102	vacuum	GeneAll® GenEx ™	M for iso	lation of	total DNA wit	hout shin
GeneAll® Exgen	eTM for is	olation o	f total DNA		General Gener	101 1301			nout spin c
		100	104-101	spin /	GenEx [™] Blood	Sx	500	220-101	solution
	mini	250	104-152	vacuum	Geriez Blood	Lx	100	220-103	solution
T 0.4		26	104-226	spin /			100	221-101	301411011
Tissue SV	Midi	Midi 100	104-201	vacuum	GenEx [™] Cell	Sx	500	221-105	solution
	MAXI	10	104-310	spin /	-	Lx	100	221-301	solution
	MAXI	26	104-326	vacuum			100	222-101	1.7
	mini	100	109-101	spin /	$GenEx^{TM}$ Tissue	Sx	500	222-105	solution
		250	109-152	vacuum	_	Lx	100	222-301	solution
Tissue plus! SV	Midi	26	109-226	spin /					
F :		100	109-201	vacuum					
				spin /					
	MAXI	10	109-310	spin / vacuum					

Ordering Information

Produ	cts	Scale	Size	Cat. No.	lype				
GeneAll® GenEx TM for isolation of total DNA									
		Sx	100	227-101					
GenEx [™] Pla	ant	Mx	100	227-201	solution				
-		Lx	100	227-301					
		Sx	100	228-101					
GenEx [™] Pla	ant plus!	Mx	50	228-250	solution				
		Lx	20	228-320					

Duadwate Carla Cias C (N

GeneAll® *DirEx™* series

 $\mathsf{Dir}\mathsf{Ex}^{\mathsf{TM}}\mathit{Fast}\mathsf{-}\mathsf{Buccal}\;\mathsf{swab}$

DirEx[™] Fast-Cigarette

	'	,	
DirEx [™]	100	250-101	solution
DirEx [™] Fast-Tissue	96 T	260-011	solution
DirEx [™] Fast-Cultured cell	96 T	260-021	solution
DirEx [™] Fast-Whole blood	96 T	260-031	solution
DirEx [™] Fast-Blood stain	96 T	260-041	solution
DirEx [™] Fast-Hair	96 T	260-051	solution

96 T

96 T

260-061

260-071

solution

solution

for preparation of PCR-template without extraction

GeneAll® RNA series for preperation of total RNA

RiboEx [™]	mini	100	301-001	solution
NIDUEX	mini	200	301-002	SOIUIION
Hybrid-R [™]	mini	100	305-101	spin
Hybrid-R [™] Blood RNA	\mini	50	315-150	spin
$\overline{\text{Hybrid-R}^{\text{TM}} \text{miRNA}}$	mini	50	325-150	spin
RiboEx [™] LS	mini	100	302-001	solution
KIDOEX LS	mini	200	302-002	SOIUUON
$Riboclear^{TM}$	mini	50	303-150	spin
Riboclear [™] plus!	mini	50	313-150	spin
Ribospin™	mini	50	304-150	spin
Du · TMu	mini	50	314-150	onin
Ribospin [™] II		300	314-103	spin
Ribospin [™] vRD	mini	50	302-150	spin
Ribospin [™] vRD <i>plus!</i>	mini	50	312-150	spin
Ribospin [™] vRD II	mini	50	322-150	spin
Ribospin [™] Plant	mini	50	307-150	spin
Ribospin [™] Seed / Fruit	mini	50	317-150	spin
Allspin TM	mini	50	306-150	spin
RiboSaver™	mini	100	351-001	solution

	Scale	SIZE	Cat. 140.	Type
GeneAll® AmpO	NF TM for	PCR an	nhlification	
selleAll 7m.po	10.	250 U		
Tag DNA polymeras	e	500 U	501-050	(2.5 U/ µl)
, , ,		I,000 U	501-100	- (
		250 U	502-025	
lpha-Taq DNA polyme	rase	500 U	502-050	(2.5 ∪/µℓ)
		I,000 U	502-100	
		250 U	504-025	
lpha-Pfu DNA polyme	rase	500 U	504-050	(2.5 U/µℓ)
		I,000 U	504-100	
		250 U	505-025	
Fast-Pfu DNA polymerase		500 U	505-050	(2.5 U/µℓ)
polymerase		1,000 U	505-100	
		250 U	531-025	
Hotstart Taq DNA polymerase		500 U	531-050	(2.5 U/µℓ)
		1,000 U	531-100	
		20 µl	521-200	1 122 1
T D :		50 µl	521-500	- lyophilized
Taq Premix	96 tubes	20 µl	526-200	1.0
		50 µl	526-500	- solution
		20 µl	522-200	1 12 1
or Tra Describe	0/+	50 µl	522-500	- lyophilized
lpha-Taq Premix	96 tubes	20 µl	527-200	
		50 µl	527-500	- solution
		20 µl	525-200	1.6
HS-Taq Premix	96 tubes	50 μl	525-500	solution
		20 µl	520-200	lyophilized
	96 tubes	50 μl	523-500	solution
Taq Premix (w/o dye)	96 tubes	s 20 µl	524-200	lyophilized
dNTPs mix		500 µl	509-020	2.5 mM eac
dNTPs set (set of dATP, dCTP, dGTP and	d dTTP)	I ml x 4 tubes	509-040	100 mM

Products Scale Size Cat. No. Type

^{*} Each dNTPs is available

Ordering Information

Products	Scale	Size	Cat. No.	Туре	Products	Size	Cat. No
GeneAll® A mp/	Master™	for PCR	amplification		GeneAll® Protein series		
	0.5 ml x 2	tubes	541-010	solution	ProtinEx TM	701-001	solution
Taq Master mix	0.5 ml x l	0 tubes	541-050	solution	Animal cell / tissue	701-001	SOIULIOI
Of Transferrance	0.5 ml x 2	tubes	542-010	solution	PAGESTA™		
α -Taq Master mix	0.5 ml x l	0 tubes	542-050	solution	Reducing I ml × 10 tubes	751-001	solution
LICT- Master as	0.5 ml x 2	tubes	545-010	solution	Sample Buffer		
HS-Taq Master mix	0.5 ml x l	0 tubes	545-050	solution	• TM		
DC N4	0.5 ml x 2	tubes	543-010	solution	GeneAll® STEADi™ for auton	natic nucleic o	ncid puriticat
lpha-Pfu Master mix	0.5 ml x l	0 tubes	543-050	solution	STEADi [™] 12 Instrument		GST012
- AII®	• . TM				STEADi [™] 24 Instrument		GST024
GeneAll® Hype					STEADi™ Genomic DNA	96	401-104
Reverse Transcripta			601-100	solution	Cell / Tissue Kit	70	T01-10T
RT Master mix	0.5 ml ×	2 tubes	601-710	solution	STEADi [™] Genomic DNA Blood K	it 96	402-105
RT Master mix with oligo (dT) ₂₀	0.5 ml ×	2 tubes	601-730	solution	STEADi™ Bacteria DNA Kit	96	403-106
RT Master mix with random hexamer	0.5 ml ×	2 tubes	601-740	solution	STEADi [™] Total RNA Kit STEADi [™] Viral DNA / RNA Kit	96 96	404-304
RT Premix	96 tubes,	20 µl	601-602	solution			
RT Premix with oligo (dT) ₂₀	96 tubes,	20 µl	601-632	solution	STEADi [™] CFC Seed DNA / RNA k	it 96	406-C02
RT Premix with random hexamer	96 tubes,	20 µl	601-642	solution			
One-step RT-PCR Master mix	0.5 ml ×	2 tubes	602-110	solution			
One-step RT-PCR Premix	96 tubes,	20 µl	602-102	solution			
First strand Synthesis Kit	50 rea	action	605-005	solution			
ZymAll [™] RNase Inhibitor	10,00	00 U	605-010	solution			
ZymAll [™] RNase Inhibitor	4,00)0 U	605-004	solution			
GeneAll® Real	Amp TM for	aPCR a	mplification				
SYBR qPCR Master		20 μl	801-020				
mix (2X, Low ROX)		20 με	801-050	solution			
SYBR qPCR Master		20 μl	801-021				
mix (2X, High ROX)		20 μl	801-051	solution			



www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro, Songpa-gu, Seoul, South Korea 138-859

E-mail : sales@geneall.com

Tel: 82-2-407-0096 Fax: 82-2-407-0779

©2016 GeneAll Biotechnology, All right reserved

Edited by **KEJ**Designed by **KDY**2016.07