Ribospin[™] vRD (Cat. No. 302-150, 302-103) / Ribospin[™] vRD II (Cat. No. 322-150, 322-103)

Purification of viral DNA/RNA using Ribospin[™] vRD and Ribospin[™] vRD II

These protocols are especially designed for isolation of viral DNA/RNA from whole blood, saliva, urine, stool, tissue, sputum and dried blood spot

Preparation

- Preheat a water bath at 56℃, 85℃
- Proteinase K (20 mg/ml)
- 1.5 ml microcentrifuge tube
- 1X PBS
- 2 ml Glass Bead or 5 mm diameter Stainless Steel Bead
- Buffer CL

Whole Blood, Urine

- 1. Dilute whole blood to 1:5 volume ratio with 1X PBS.
- 2. Continue with step 1 of each protocol.

2 Stool

* Good performance can be achieved without pre-treatment, but pre-treatment described below might be an option to enhance extraction efficiency.

- Add 50~100 mg of stool or 50~100 μl of diarrhea stool to 2 ml Glass Bead tube or 5 mm diameter Stainless Steel Bead tube.
- Add 400 µl of Buffer VL to the tube. (Ribospin[™] vRD II : Buffer NVL) Buffer VL of Ribospin[™] vRD or Buffer NVL of Ribospin vRD[™] II can be replaced with PBS.
- 3. Vortex for more than 1 min.
- 4. Incubate for 1 min at room temperature.
- 5. Transfer 200 μl of supernatant to new tube without disturbing the precipitate.
- 6. Continue with step 1 of each protocol.

3 Tissue

* TissueLyser II method

- 1. Add 50~100 mg of tissue sample to 2 ml tube with 6 mm diameter Stainless Steel Bead.
- 2. Add 400 µl of Buffer VL to the tube. (Ribospin[™] vRD II : Buffer NVL)
- 3. Grind the tissue at 30 Hz for 30 sec using TissueLyser II.
- 4. Add 400 µl of Buffer VL. (Ribospin[™] vRD II : Buffer NVL)
- 5. Vortex for 30 sec and incubate for 1 min at room temperature.
- 6. Transfer 200 μl of supernatant to new tube without disturbing the precipitate.
- 7. Continue with step 1 of each protocol.

Supplementary Protocol

Ribospin[™] vRD (Cat. No. 302-150, 302-103) / Ribospin[™] vRD II (Cat. No. 322-150, 322-103)

* Pestle method

- 1. Grind 100 mg of tissue with Pestle for less than 1 min.
- 2. Add 400 µl of Buffer VL. (Ribospin[™] vRD II : Buffer NVL)
- 3. Vortex for 30 sec and centrifuge 13,000 rpm for 1 min at room temperature.
- 4. Transfer 200 μl of supernatant to new tube without disturbing the precipitate.
- 5. Continue with step 1 of each protocol.

4 Sputum, Saliva

- 1. Dilute the sputum to 1:1 volume ratio with PBS.
- 2. Vortex for 1 or 5 min depending on the sputum viscosity.
- 3. Continue with step 1 of each protocol.

5 Dried blood spot (FTA card)

- Place 1~3 punched-out circles from a dried blood spot into a 1.5 ml microcentrifuge tube and add 300 µl of Buffer VL. (Ribospin[™] vRD II : Buffer NVL)
- Incubate at 85°C for 10 min.
 (Optional) Add 20 µl of Proteinase K solution (20 mg/ml, not provided) and 20 µl of 1M DTT, vortex to mix, and incubate at 56°C for 10 min.
- 3. Centrifuge 13,000 rpm for 1 min at room temperature.
- 4. Transfer 200 µl of supernatant to new tube without disturbing the precipitate.
- 5. Add 300 μl of Buffer RB1 and gently mix or inverting. Incubate at room temperature for 1 min.
- 6. Transfer all of the mixture, 500 µl to a Column Type V (mini). [Ribospin[™] vRD II : Column Type S (micro)]
- 7. Continue with step 6 of Ribospin[™] vRD protocol. (Ribospin[™] vRD II : continue with step 7)

Purification of viral DNA/RNA from whole blood, saliva, urine, stool, tissue, sputum and dried blood spot using Ribospin[™] vRD and Ribospin[™] vRD II.