# Allspin<sup>™</sup>

For total DNA & RNA isolation from tissues and cultured cells

TOTAL DNA/RNA PURIFICATION HANDBOOK



## **Customer & Technical Support**

Should you have any further questions, do not hesitate to contact us. We appreciate your comments and advice.

## **Contact Information**

www.geneall.com Tel : 82-2-407-0096 Fax : 82-2-407-0779 E-mail(Order/Sales) : sales@geneall.com E-mail(Tech. Info.) : tech@geneall.com

# Visit GeneAll<sup>®</sup> Community

www.geneall.com www.geneall.co.kr

## Trademarks

Hybrid-Q<sup>™</sup>, Exprep<sup>™</sup>, Exfection<sup>™</sup>, Expin<sup>™</sup>, Exgene<sup>™</sup>, *GenEx<sup>™</sup>*, DirEx<sup>™</sup>, Hybrid-R<sup>™</sup>, RiboEx<sup>™</sup>, Ribospin<sup>™</sup>, Riboclear<sup>™</sup>, Allspin<sup>™</sup>, Allspin<sup>™</sup>, RiboSaver<sup>™</sup>, EzClear<sup>™</sup>, EzSep<sup>™</sup>, EzPure<sup>™</sup>, EzPass<sup>™</sup>, AmpONE<sup>™</sup>, AmpMaster<sup>™</sup>, RealAmp<sup>™</sup>, HyperScript<sup>™</sup>, ProtinEx<sup>™</sup>, PAGESTA<sup>™</sup>, STEADi<sup>™</sup> are trademarks of GeneAll Biotechnology Co., Ltd.


© 2017 GeneAll Biotechnology, all rights reserved.


This protocol handbook is included in :

GeneAll<sup>®</sup> Allspin<sup>™</sup> (306-150)

Visit www.geneall.com or www.geneall.co.kr for FAQ, Q&A and more information.

# **Brief Protocol**





# INDEX

Brief Protocol	03
Index	05

- Kit Contents 06
- Materials Not Provided
- Product Specifications 07 Quality Control Storage Conditions
- Safety Information 08 Product Disclaimer

Prevention of RNase Contamination

- Sample Amount and Average Yield 09
  - Product Description 10

### Protocol

- Cultured animal cells
  - Animal tissues 15
- Troubleshooting Guide 20
  - Appendix I 23
  - Appendix 2 27
  - Appendix 3 28

### Ordering Information 30

# **Kit Contents**

Cat. No.	306-150	Storage
Components	Quantity	Slorage
Buffer CTL	45 ml	
Buffer GW1	45 ml	
Buffer BW (concentrate) *	16 ml	
Buffer RNW (concentrate) * †	12 ml x 2	Room
Buffer AE **	15 ml	temperature
Nuclease-free water	15 ml	(15~25°C)
Mini column type B (with collection tube)	50	
Column type W (mini) (with collection tube)	50	
2.0 ml collection tube	50	
1.5 ml microcentrifuge tube	100	

\* Before first use, add absolute ethanol (ACS grade or better) into Buffer BW and RNW as indicated on the bottle.

+ Contains sodium azide as a preservative

\*\* 10 mM TrisCl, pH 9.0, 0.5 mM EDTA

## **Materials Not Provided**

#### Reagents

- Tissue storage buffer to protect RNA from RNase
- Absolute ethanol, ACS grade or better

#### Disposable materials

- RNase-free pipet tips
- Disposable gloves

### Equipments

- Equipment for disrupting sample
- Microcentrifuge, Vortex mixer
- Suitable protector (ex; lab coat, goggles, etc)

Allspin™ total DNA / RNA purification kit						
Specification	Column type B for DNA	Column type W for RNA				
Туре	Spin / mini	Spin / mini				
Maximum amount of starting samples	30 mg / prep or ~1 x 10 <sup>7</sup> cells / prep	30 mg / prep or $\sim 1 \times 10^7$ cells / prep				
Preparation time	≥30 min	≥30 min				
Maximum loading volume of mini column	750 <i>µ</i> I	750 <i>µ</i> I				
Minimum elution volume	50 <i>µ</i> I	30 <i>µ</i> I				

## **Quality Control**

All components in GeneAll<sup>®</sup> Allspin<sup>™</sup> are manufactured in strictly clean conditions, and its degree of cleanness is monitored periodically. Quality control is carried out thoroughly from lot to lot, and only the qualified kits are approved to be delivered.

## **Storage Conditions**

All components of GeneAll<sup>®</sup> Allspin<sup>TM</sup> should be stored at room temperature ( $15\sim25^{\circ}$ C). It should be protected from exposure to direct sunlight.

During shipment or storage under cool ambient condition, a precipitate can be formed in Buffer CTL, GWI, and BW. In such a case, heat the bottle to  $50^{\circ}$ C to dissolve completely. Using precipitated buffers will lead to poor DNA recovery. GeneAll<sup>®</sup> Allspin<sup>TM</sup> is guaranteed until the expiration date printed on the product box.

## **Safety Information**

The buffers included in the Allspin<sup>TM</sup> total DNA / RNA purificaion kit contain irritants which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken when handling such materials. Always wear gloves and eye protection, and follow standard safety precautions.

Buffer CTL, GWI, and BW contain chaotropic agents, which can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

## **Product Disclaimer**

GeneAll<sup>®</sup> Allspin<sup>™</sup> is for research use only, not for use in diagnostic procedure.

## **Preventing RNase Contamination**

RNase can be introduced accidentally during RNA purification. Wear disposable gloves always, because skin often contains bacteria and molds that can be a source of RNase contamination. Use sterile, disposable plastic wares and automatic pipettes to prevent cross-contamination of RNase from shared equipment.

## Sample Amount and Average Yield

Amount of starting material	Sample type	Average yield of genomic DNA	Average yield of total RNA
Cultured cells	СНО	$\sim$ 7 $\mu$ g	~15 µg
$(=   \times   0^6)$	RAW264.7	~10 µg	~20 µg
Tissue (rat)	Liver	$\sim$ 25 $\mu$ g	~60 µg
(10 mg / prep)	Kidney	$\sim$ 25 $\mu$ g	~30 µg
	Brain	~ 2 µg	~10 µg
	Heart	~10 µg	~9 µg
	Spleen	~70 µg	$\sim$ 80 $\mu$ g

The yield of genomic DNA and total RNA may vary depending on the type of tissue or cells from which it is obtained.

Allspin<sup>™</sup> total DNA / RNA purification kit provides a convenient method for the isolation of total DNA and total RNA simultaneously from a single sample of tissue or cultured cells. It utilizes the optimized buffer system and the advanced silica-binding technology to purify nucleic acid sufficiently pure for many applications, instead of conventional alcohol precipitation and phenol / chloroform extraction.

DNA and RNA are purified separately from a same sample but successive procedure using Column Type B and W respectively. Alternatively, both DNA and RNA can be copurified into a single tube by the modified procedure at appendix 1. Whole procedure can be performed in just 30 minutes and the length of obtained DNA is up to 50 kb (average is 30 kb) and that of RNA is longer than 200 nucleotides.

To obtain pure DNA, samples are homogenized in Buffer CTL, containing guanidine salt, which rapidly lyses cells and inactivates nucleases and then the lysate is applied into Column type B. During centrifugation, DNA is bound specifically to column membrane and RNA is passed through it and goes into the collection tube. The membrane is washed by a series of washing buffers and bound DNA is eluted by Buffer AE.

For purification of pure RNA, the passed-through which is obtained from fraction of DNA purification procedure is mixed with ethanol and this mixture is loaded into column W. RNA is bound specifically during centrifugation. After washed by Buffer GWI and RNW, pure RNA is eluted by Nuclease-free water.

Purified DNA and RNA are fully suitable for the isolation of Poly  $A^+$  RNA, southern and northern blotting, dot blotting, in vitro transcription, cloning, RT-PCR and other analytical procedures.

# Allspin<sup>TM</sup> total DNA / RNA purificaion kit

## PROTOCOL for cultured animal cells

### Before experiment

- \* Before first use, add absolute ethanol (ACS grade or better) to Buffer BW, Buffer RNW as indicated on the bottle.
- \* Prepare 70% ethanol
- \* All centrifugation should be carried out at room temperature in a microcentrifuge.

### I. Harvest cell samples in a 1.5 ml microcentrifuge tube.

### Cells grown in monolayer

Harvest  $5 \times 10^6$  cells carefully using scraper. And pellet cells by centrifugation at low speed (below 3000 rpm) for 5 min, then discard the culture medium.

### Cells grown in suspension

Pellet  $5 \times 10^6$  cells by centrifugation at low speed (below 3000 rpm) for 5 min, then discard the culture medium.

\* Do not wash cells before lysing with Buffer CTL as this may contribute to mRNA degradation.

2. Add 350  $\mu$ I of Buffer CTL to the tube and lyse the sample by pipetting or homogenizing.

Lyse the 5 x 10<sup>6</sup> cells in 350  $\mu$ I Buffer CTL using pipetting. An insufficient lysis may result in low RNA recovery rate or mini column clogging.

3. Incubate the lysate for 3 min at room temperature.

This step allows nucleoprotein complexes to completely dissociate.

# Cell



**12** GeneAll<sup>®</sup> Allspin<sup>™</sup> Protocol Handbook

- Transfer the lysate to a Column Type B (red ring). Centrifuge at ≥ 10,000 x g for 30 sec at room temperature. Maximun volume is 750 µl and DNA is bound to membrane through this step.
- 5. Transfer the mini column to a new 2.0 ml collection tube (provided), and store at room temperature. Use the pass-through for total RNA purification.

Make sure that no lysate remains in the mini column after centrifugation. If the residual lysate has remained, centrifuge again at higher speed until all of the solution has passed through.

Go on to step 6 for total RNA purification from cell samples.

Go on to step 17 for genomic DNA purification from cell samples.

Total RNA purification (Blue ring column step)

Cell

- 6. Add I volume (usually 350 μl) of 70% ethanol (not provided) to the collection tube including the pass-through (step 5), and mix well by pipetting.
  Do not centrifuge.
- 7. Transfer the mixture (approximately 700  $\mu$ l) including any precipitate to a Column Type W (blue ring).
- 8. Centrifuge at  $\geq 10,000 \times g$  for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.







- 9. Add 700  $\mu$ l of Buffer GWI to the mini column.
- 10. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- II. Add 500  $\mu$ I of Buffer RNW to the mini column.
- 12. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 13. Repeat the step 11 and 12.
- 14. Centrifuge at  $\geq 10,000 \text{ x}$  g for an additional 1 min at room temperature to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of Buffer RNW.

15. Add 50  $\mu$ l of Nuclease-free water to the center of the membrane in the mini column.

According to the expected yield, the volume of eluent can be adjusted.

16. Centrifuge at ≥ 10,000 x g for I min at room temperature. Purified RNA can be stored at 4°C for immediate analysis and can be stored at -70°C for long term storage.



Nuclease-free water

Ready for use!

Cel



- 17. Add 500  $\mu$ l of Buffer BW to the mini column.
- 18. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 19. Add 500  $\mu$ l of Buffer RNW to the mini column.
- 20. Centrifuge at  $\geq 10,000 \text{ x}$  g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
  - 21. Centrifuge at  $\geq$  10,000 x g for an additional 1 min at room temperature to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of Buffer RNW.

22. Add 50  $\mu$ I of Buffer AE to the center of the membrane in the mini column. Let it stand for I min.

According to the expected yield, the volume of eluent can be adjusted.

23. Centrifuge at  $\geq$  10,000 x g for 1 min at room temperature.

For long-term storage, eluting in Buffer AE is recommended. But, EDTA included in Buffer AE can inhibit some delicate enzymatic reaction, so you can avoid such latent problems by using distilled water (>pH 7.0) or Tris-Cl (>pH 8.5). When using water for elution, make sure the pH of water is higher than 7.0.



# Allspin<sup>TM</sup> total DNA / RNA purificaion kit

# PROTOCOL for animal tissues

#### Before experiment

- \* Before first use, add absolute ethanol (ACS grade or better) to Buffer BW, Buffer RNW as indicated on the bottle.
- \* Prepare 70% ethanol
- \* All centrifugation should be carried out at room temperature in a microcentrifuge.
- I. Harvest tissue samples in a 1.5 ml microcentrifuge tube. Harvest ~20 mg tissue samples in a 1.5 ml microcentrifuge tube. The recommended method for sample handling is to put directly removed fresh tissue into tissue storage buffer or to freeze the tissue rapidly in liquid nitrogen.
- 2. Add 350  $\mu$ l of Buffer CTL to the tube and disrupt and homogenize the sample by homogenization.

Before starting, add ß-mercaptoethanol (10  $\mu I$  per 1 ml) to Buffer CTL.

Homogenize  $\sim$ 20 mg of tissue samples in 350  $\mu$ l Buffer CTL using homogenizer (rotor-stator homogenizer, mortar and pestle, or bead-beater). Thoroughly disrupt the tissue in Buffer CTL and lyse the samples completely. Not clarified sample may cause clogging of the mini column in subsequent steps.

### 3. Incubate the lysate for 10 min at room temperature.

This step allows nucleoprotein complexes to completely dissociate and reduces the generated foam.

### lissue

Buffer CTL

350,µl

ZZZZZZ

0 min

4. Centrifuge at maximum speed for 3 min at room temperature and carefully transfer the supernatant to a Column Type B (red ring).

This step can help avoid clogging of the mini column caused by not clarified insoluble particles.

Maximun volume is 700  $\mu$ l and DNA is bound to membrane through this step.

 Centrifuge at ≥ 10,000 x g for 30 sec at room temperature. Transfer the mini column to a new 2.0 ml collection tube (provided), and store at room temperature. Use the passthrough for total RNA purification.

Make sure that no lysate remains in the mini column after centrifugation. If the residual lysate has remained, centrifuge again at higher speed until all of the solution has passed through.

Go on to step 6 for total RNA purification from tissue sample.

Go on to step 17 for genomic DNA purification from tissue sample.

Total RNA purification (Blue ring column step)

- 6. Add I volume (usually 350  $\mu$ I) of 70% ethanol to the collection tube including the pass-through (step 5), and mix well by pipetting. Do not centrifuge.
- 7. Transfer the mixture (approximately 700  $\mu$ l) including any precipitate to a Column Type W (blue ring).



Supernatant







- 8. Centrifuge at  $\geq 10,000 \times g$  for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 9. Add 700  $\mu$ I of Buffer GWI to the mini column.
- 10. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- II. Add 500  $\mu$ I of Buffer RNW to the mini column.
- 12. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 13. Repeat the step 11 and 12.
- 14. Centrifuge at  $\geq 10,000 \text{ x}$  g for an additional 1 min at room temperature to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of Buffer RNW.

15. Add 50  $\mu$ l of Nuclease-free water to the center of the membrane in the mini column.

According to the expected yield, the volume of eluent can be adjusted.





16. Centrifuge at ≥ 10,000 x g for 1 min at room temperature. Purified RNA can be stored at 4°C for immediate analysis and can be stored at -70°C for long term storage.

In tissue samples, occasionally genomic DNA contamination is happened in the RNA eluate because of lots of sample amount. In this case, reduce the starting sample amounts or use a DNase I following the manufacture's instruction manual.

Genomic DNA purification (Red ring column step)

- **Tissue** 17. Add 500  $\mu$ l of Buffer BW to the mini column.
  - 18. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
  - 19. Add 500  $\mu$ l of Buffer RNW to the mini column.
  - 20. Centrifuge at  $\geq 10,000 \times g$  for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
  - 21. Centrifuge at  $\geq$  10,000 x g for an additional 1 min at room temperature to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of Buffer RNW.








22. Add 100  $\mu$ l of Buffer AE to the center of the membrane in the mini column. Let it stand for I min.

According to the expected yield, the volume of eluent can be adjusted.

23. Centrifuge at  $\geq$  10,000 x g for 1 min at room temperature.

For long-term storage, eluting in Buffer AE is recommended. But, EDTA included in Buffer AE can inhibit some delicate enzymatic reaction, so you can avoid such latent problems by using distilled water (>pH 7.0) or Tris-Cl (>pH 8.5). When using water for elution, make sure the pH of water is higher than 7.0.



# **Troubleshooting Guide**

Facts	Possible Causes	Suggestions		
Low yield	Insufficient homogenization of the sample with Buffer CTL	Be sure to incubate after homogenization. Confirm the completely homogenized sample in Buffer CTL.		
	Too much starting sample	Reduce the amount of starting sample. Especially for tissue sample, use the correct amount of starting sample.		
	Poor quality of starting material	<b>g</b> Process the sample immediately after harves from animal. To process later, freeze the tissue rapidly in liquid nitrogen.		
	Residual culture media in the sample	Remaining culture media affect lysis and binding condition of Buffer CTL.		
Column clogging	Insufficient homogenization of the sample with Buffer CTL	Be sure to incubate after homogenization. Confirm the completely lysed sample in Buffer CTL.		
	Too much starting sample	Reduce the amount of starting sample. Especially for tissue sample, use the correct amount of starting sample.		
Low A <sub>260</sub> / A <sub>280</sub>	Insufficient homogenization of the sample with Buffer CTL	Be sure to incubate after homogenization. Confirm the completely lysed sample in Buffer CTL.		
High A260 / A280 in DNA eluate	RNA contamination	RNA may inhibit some downstream enzymatic reactions. If RNA-free DNA is required, RNase treatment should be performed following the manufacture's protocol.		
		Use Buffer AE to dilute extracted nucleic acid and to make blank before measure purity, If Buffer AE has been used for elution.		

# Troubleshooting Guide

Facts	Possible Causes	Suggestions			
DNA contamination in RNA eluate	Too much starting sample	Reduce the amount of starting sample. Especially for tissue sample, use the correct amount of starting sample.			
	No treatment Buffer GWI	Follow the Allspin <sup>™</sup> total DNA / RNA purification kit manual.			
	High DNA mass in sample	Especially, some tissue samples have high DNA content. In this case, some DNA will pass the Column Type B. And then bind to the next RNA binding step at Column Type W with RNA. Reduce the amount of starting sample or perform DNase digestion at the RNA eluate.			
RNA contamination in DNA eluate	Too much manipulated sample before process	Reduce the amount of starting sample. Especially for tissue sample, use the correct amount of starting sample.			
RNA degradation	Too much manipulated sample before process	Process the sample immediately after harvest from animal. For cultured cell, minimize washing steps in cell harvest.			
	Improper storage of RNA	Store isolated RNA at -70°C, Do not store at -20°C.			
	Use of not RNase-free reagent or disposable products	Make sure to use RNase-free products only.			
DNA degradation	Too old or mis-stored starting sample	Too old or mis-stored sample often yield degraded DNA. Use fresh sample.			
	Vigorous homogenization	Vigorous handling after addition of Buffer CTL can lead to irreversible denaturation of genomic DNA. Minimize the homogenization time and gently homogenize the sample.			

## **Troubleshooting Guide**

Facts	Possible Causes	Suggestions
Eluate does not perform well in downstream application	Residual ethanol in eluate	To remove any residual ethanol included in Buffer RNW from mini column membrane, centrifuge again for complete removal of ethanol.
	Use of Buffer BW and RNW in the wrong order	Ensure that Buffer BW and RNW are used in the correct order in the protocol. If used in the wrong order, perform the last washing step with Buffer RNW.

# **APPENDIX** . Copurification of total DNA and RNA in a single tube

Allspin<sup>™</sup> total DNA / RNA purificaion kit is convenient kit for isolation of genomic DNA and total RNA from tissues or cultured cells. Allspin<sup>™</sup> total DNA / RNA purificaion kit is designed to purify DNA and RNA separately, but, DNA and RNA can be purified simultaneously in a single tube using the modified method. After lysis, the lysate is mixed with ethanol then apply to Column Type B to bind DNA and RNA on the membrane. And then the impurities on the membrane are washed away by two different wash buffers. At last pure RNA and DNA are eluted by Nuclease-free water. The eluate should be treated with care because RNA is very sensitive to contamination.

### Protocol for cultured cell samples.

I. Harvest cell samples in a 1.5 ml microcentrifuge tube.

Cells grown in monolayer

Harvest  $5 \times 10^6$  cells carefully using scraper. And pellet cells by centrifugation at low speed (below 3000 rpm) for 5 min, then discard the culture medium.

### Cells grown in suspension

Pellet 5 x 10 $^{6}$  cells by centrifugation at low speed (below 3000 rpm) for 5 min, then discard the culture medium.

\* Do not wash cells before lysing with Buffer CTL as this may contribute to mRNA degradation.

# 2. Add 350 $\mu$ l of Buffer CTL to the tube and lyse the sample by pipetting or homogenizing.

Lyse the 5 x 10° cells in 350  $\mu$ l Buffer CTL using pipetting. An insufficient lysis may result in low RNA recovery rate or column clogging.

### 3. Incubate the lysate for 3 min at room temperature.

This step allows nucleoprotein complexes to completely dissociate.

- 4. Add I volume (usually  $350 \mu$ I) of absolute ethanol to the lysate, and mix well by pipetting or vortexing. Do not centrifuge.
- 5. Transfer the mixture (approximately 700  $\mu$ I) including any precipitate to a Column Type B (red ring).
- Centrifuge at ≥ 10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 7. Add 500  $\mu$ I of Buffer BW to the mini column.
- Centrifuge at ≥10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 9. Add 500 µl of Buffer RNW to the mini column.
- 10. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 11. Centrifuge at  $\geq$  10,000 x g for an additional 1 min at room temperature to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of Buffer RNW.

12. Add 100  $\mu$ l of Nuclease-free water to the center of the membrane in the mini column. Let it stand for I min.

According to the expected yield, the volume of eluent can be adjusted.

13. Centrifuge at  $\geq$  10,000 x g for 1 min at room temperature.

### Protocol for tissue samples.

### I. Harvest tissue samples in a 1.5 ml microcentrifuge tube.

Harvest  $\sim$ 20 mg tissue samples in a 1.5 ml microcentrifuge tube. The recommended method for sample handling is to put directly removed fresh tissue into tissue storage buffer or to freeze the tissue rapidly in liquid nitrogen.

# 2. Add 350 $\mu$ l of Buffer CTL to the tube and disrupt and homogenize the sample by homogenization.

Before starting, add ß-mercaptoethanol (10  $\mu$ l per 1 ml) to Buffer CTL.

Homogenize  $\sim 20$  mg of tissue samples in 350  $\mu$ l Buffer CTL using homogenizer (rotor-stator homogenizer, mortar and pestle, or bead-beater). Thoroughly disrupt the tissue in Buffer CTL and lyse the samples completely. Not clarified sample may cause clogging of the mini column in subsequent steps.

### 3. Incubate the lysate for 10 min at room temperature.

This step allows nucleoprotein complexes to completely dissociate and reduces the generated foam.

4. Centrifuge at maximum speed for 3 min at room temperature and carefully transfer the supernatant to a new 1.5 ml microcentrifuge tube (not provided).

This step can help avoid clogging of the mini column caused by not clarified insoluble particles.

Maximun volume is 700  $\mu$ l and DNA is bound to membrane through this step.

# 5. Add I volume (usually 350 $\mu$ I) of absolute ethanol to the lysate, and mix well by pipetting or vortexing. Do not centrifuge.

- 6. Transfer the mixture (approximately 700  $\mu$ l) including any precipitate to a Column Type B (red ring).
- Centrifuge at ≥10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 8. Add 500  $\mu$ l of Buffer BW to the mini column.

- Centrifuge at ≥ 10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 10. Add 500  $\mu$ l of Buffer RNW to the mini column.
- 11. Centrifuge at  $\geq$  10,000 x g for 30 sec at room temperature. Discard the pass-through and reinsert the mini column back into the collection tube.
- 12. Centrifuge at  $\geq$  10,000 x g for an additional I min at room temperature to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of Buffer RNW.

13. Add 100  $\mu$ l of Nuclease-free water to the center of the membrane in the mini column. Let it stand for I min.

According to the expected yield, the volume of eluent can be adjusted.

14. Centrifuge at  $\geq$  10,000 x g for 1 min at room temperature.

# $\begin{array}{c} \text{APPENDIX 2.} \\ \text{Confirmation of RNA yield and purity by} \\ \text{UV absorbance} \end{array}$

### **Concentration of RNA**

The concentration of RNA can be determined by the absorbance at 260 nm using spectrophotometer. For the convenient measurement, we recommend using the NanoDrop<sup>®</sup> which can reduce your RNA sample and time. If not, you need to dilute the RNA samples to measure the concentration through traditional spectrophotometer. The value of  $A_{260}$  should be between 0.15 and 1.00. Be sure to calibrate the spectrophotometer with the same solution used for dilution.

An absorbance of 1 at 260 nm is correspond to about 40  $\mu$ g RNA / ml at a neutral pH. Therefore, the concentration of RNA was calculated by the formula shown below.

A\_{260} X dilution factor X 40 = RNA  $\mu$ g / ml

### **Purity of RNA**

To confirm the RNA purity, you should read the ratio of  $A_{260}$  /  $A_{280}.$  Pure RNA is in the range of  $1.8{\sim}2.2.$ 

### APPENDIX 3. Formaldehyde agarose gel electrophoresis (Denaturing gel method)

A denaturing agarose gel is routinely used for the assessment of the quality of an RNA preparation. After preparation, RNA forms secondary structure via intramolecular base pairing. Therefore, it is very difficult to get the exact result of electrophoresis because of migrating inaccuracy. However, the denaturing gel denatures the secondary structure of RNA and makes an accurate migration. To confirm the RNA band, the gel should be transferred to a UV transilluminator after electrophoresis. Mainly, two RNA bands are shown. In case of animal sample, the 28S and 18S rRNA bands are confirmed on the gel. If they are intact, the RNA bands should be sharp and the intensity of upper band should be about twice that of the lower band.

### Prepare the denaturing gel

- I. Put I g agarose in 72 ml water and heat to dissolve thoroughly.
- 2. Cool to 60°C.
- 3. Add 10 ml of 10X MOPS buffer, 18 ml of 37% formaldehyde, and 1  $\mu$ l of a 10 mg / ml ethidium bromide (EtBr).
- 4. Mix well then pour the gel into the gel tray and cool to solidify it.
- 5. Transfer the solidified gel from tray to tank, and add enough IX MOPS running buffer to cover the gel.

### Prepare the RNA sample

- 1. Make the mixture. ?  $\mu$ I RNA (up to 20  $\mu$ g) 2  $\mu$ I 10X MOPS electrophoresis buffer 4  $\mu$ I formaldehyde 10  $\mu$ I formamide
- 2. Incubate the mixture for 15 minutes at  $65^{\circ}$ C.
- 3. Chill the sample for 5 minutes in ice.
- 4. Add  $2 \mu l$  of 10X formaldehyde gel-loading dye to the mixture.
- 5. Load the mixture in a denaturing gel which is covered with a sufficient IX MOPS electrophoresis buffer.
- 6. Run the gel and confirm the RNA band on transilluminator.

Occasionally, gel destaining may be needed to increase the visibility of the bands of RNA in  $dH_2O$  for several hours.

### **Composition of buffers**

### - I0X MOPS buffer

0.2 M MOPS 20 mM sodium acetate 10 mM EDTA pH to 7.0 with NaOH

### - 10X formaldehyde gel-loading dye

50% glycerol 10 mM EDTA 0.25% (w / v) bromophenol blue 0.25% (w / v) xylene cyanol FF

### \* Caution

When working with these chemicals, always use gloves and eye protector to avoid contact with skin and cloth. Especially, formaldehyde and ethidium bromide (EtBr) should be handled in a fume hood.

# **Ordering Information**

Products	Scale	Size	Cat. No.	Туре	Products	Scale	Size	Cat. No.	Туре
GeneAll® <b>Hybri</b> d	<b>I-Q™</b> fo	r rapid p	reparation of	plasmid DNA	GeneAll® Exgene	тм <sub>for is</sub>	olation o	f total DNA	
Plasmid Rapidprep	mini	50	100-150	onin		mini	100	105-101	spin /
	mini	200	100-102	spin	_	111111	250	105-152	vacuum
					Blood SV	Midi	26	105-226	spin /
GeneAll® <b>Expre</b> f	<b>b<sup>™</sup></b> for p	reparatio	n of plasmid l	DNA	BIOOD 3V	TIIQI	100	105-201	vacuum
		50	101-150	spin /		MAXI	10	105-310	spin /
	mini	200	101-102	vacuum			26	105-326	vacuum
Diagonaid CV/		26	101-226	spin /		mini	100	106-101	spin /
Plasmid SV	Midi	50	101-250		Cell SV -		250	106-152	vacuum
		100	101-201	vacuum	CCII 3V	MAXI	10	106-310	spin /
GeneAll® Exfect	ion <sup>TM</sup>						26	106-326	vacuum
		transfect	ion-grade pla	smid DNA		mini	100	108-101	spin /
		50	- 50	spin /	-		250	108-152	vacuum
Plasmid LE	mini	200	- 02	vacuum	Clinic SV	Midi	26	108-226	spin /
(Low Endotoxin)		26	-226	spin /	CIIIIC SV	DILLI	100	108-201	vacuum
	Midi	100	-20	vacuum		MAXI	10	108-310	spin /
Plasmid EF	N.C. P.	20	2 -220		MAN	26	108-326	vacuum	
(Endotoxin Free)	Midi	100	2 -20	spin	Genomic DNA micro	)	50	8-050	spin
							100	7- 0	spin /
GeneAll® <b>Expin</b> <sup>™</sup>	м for bur	ification	of fragment D	NA		mini	250	7- 52	vacuum
1		50	102-150	spin /	-	N.C. P.	26	7-226	spin /
Gel SV	mini	200	02-150 spin / Plant SV  02-102 vacuum	Midi	100	7-20	vacuum		
		50	103-150	spin /	-		10	7-3 0	spin /
PCR SV	mini	200	103-102	vacuum		MAXI	26	117-326	vacuum
		50	113-150	spin /	Soil DNA mini	mini	50	4- 50	spin
CleanUp SV	mini	200	113-102	vacuum	Stool DNA mini	mini	50	115-150	spin
			115 102						spin
		50	112-150	coin /	Viral DNA / RNA	mini	50	128-150	spin
Combo GP	mini	50 200	112-150	spin / vacuum			50 50		spin
		200	2- 02	spin / vacuum	Viral DNA / RNA FFPE Tissue DNA	mini mini		128-150	
		200 solation o	112-102 f total DNA	vacuum		mini	50 250	128-150 138-150 138-152	spin spin
		200 colation o	112-102 f total DNA 104-101	vacuum spin /	FFPE Tissue DNA	mini	50 250 lation of t	128-150 138-150 138-152 total DNA wit	spin spin
	<b>e<sup>™</sup> for is</b>	200 solation o 100 250	112-102 f total DNA 104-101 104-152	vacuum spin / vacuum	FFPE Tissue DNA	mini	50 250 lation of 1	128-150 138-150 138-152 total DNA wite 220-101	spin spin
	<b>e<sup>™</sup> for is</b>	200 colation o 100 250 26	112-102 f total DNA 104-101 104-152 104-226	spin / vacuum spin /	FFPE Tissue DNA	mini <b>m</b> for isol Sx	50 250 lation of 1 100 500	128-150 138-150 138-152 total DNA wite 220-101 220-105	spin spin thout spin co solution
GeneAll® Exgen	<b>e<sup>TM</sup></b> for is mini	200 solation o 100 250 26 100	112-102 f total DNA 104-101 104-152 104-226 104-201	spin / vacuum spin / vacuum	FFPE Tissue DNA	mini <b>m</b> for isol	50 250 <i>lation of t</i> 100 500 100	128-150 138-150 138-152 total DNA wite 220-101 220-105 220-301	spin spin thout spin co
GeneAll® Exgen	<b>e<sup>TM</sup></b> for is mini	200 solation of 100 250 26 100 10	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310	vacuum spin / vacuum spin / vacuum spin /	FFPE Tissue DNA <b>GeneAll® GenEx<sup>™</sup></b> <i>GenEx<sup>™</sup></i> Blood	mini <b>m</b> for isol Sx	50 250 <i>lation of t</i> 100 500 100	128-150 138-150 138-152 total DNA wite 220-101 220-105 220-301 221-101	spin spin thout spin co solution
GeneAll® Exgen	<b>e<sup>TM</sup></b> for is mini Midi	200 solation o 250 26 100 10 26	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326	spin / vacuum spin / vacuum spin / vacuum	FFPE Tissue DNA	mini <sup>M</sup> for isol Sx Lx Sx	50 250 lation of t 100 500 100 500	128-150 138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105	spin spin thout spin co solution solution solution
GeneAll® Exgen	<b>e<sup>TM</sup></b> for is mini Midi	200 colation o 100 250 26 100 26 100 26 100	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101	spin / vacuum spin / vacuum spin / vacuum spin / vacuum	FFPE Tissue DNA <b>GeneAll® GenEx<sup>™</sup></b> <i>GenEx<sup>™</sup></i> Blood	mini <i>for isol</i> Sx Lx	50 250 lation of t 100 500 100 500 100	128-150 138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301	spin spin thout spin co solution solution
GeneAll® Exgen	e <sup>TM</sup> for is mini Midi MAXI	200 solation of 250 26 100 10 26 100 26 100 250	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	spin / vacuum spin / vacuum spin / vacuum	FFPE Tissue DNA <b>GeneAll® GenEx<sup>™</sup></b> GenEx <sup>™</sup> Blood GenEx <sup>™</sup> Cell	mini <sup>M</sup> for isol Sx Lx Sx	50 250 100 500 100 100 500 100 100	128-150 138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101	spin spin thout spin co solution solution solution
GeneAll® Exgent	e <sup>TM</sup> for is mini Midi MAXI mini	200 solation o 250 26 100 10 26 100 250 26 26	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152 109-226	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	FFPE Tissue DNA <b>GeneAll® GenEx<sup>™</sup></b> <i>GenEx<sup>™</sup></i> Blood	mini for isol Sx Lx Sx Lx Sx Sx	50 250 100 500 100 100 500 100 100 500	128-150 138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin spin thout spin cc solution solution solution solution
GeneAll® Exgen	e <sup>TM</sup> for is mini Midi MAXI	200 solation of 250 26 100 10 26 100 26 100 250	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	FFPE Tissue DNA <b>GeneAll® GenEx<sup>™</sup></b> GenEx <sup>™</sup> Blood GenEx <sup>™</sup> Cell	mini for isol Sx Lx Sx Lx	50 250 100 500 100 100 500 100 100	128-150 138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101	spin spin thout spin co solution solution solution
GeneAll® Exgent	e <sup>TM</sup> for is mini Midi MAXI mini	200 solation o 250 26 100 10 26 100 250 26 26	112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152 109-226	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	FFPE Tissue DNA <b>GeneAll® GenEx<sup>™</sup></b> GenEx <sup>™</sup> Blood GenEx <sup>™</sup> Cell	mini for isol Sx Lx Sx Lx Sx Sx	50 250 100 500 100 100 500 100 100 500	128-150 138-150 138-152 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin spin thout spin co solution solution solution solution

Products	Scale	Size	Cat. No.	Туре
GeneAll® GenEx	for isol	lation of	total DNA	

	Sx	100	227-101	
GenEx <sup>™</sup> Plant	Mx	100	227-201	solution
	Lx	100	227-301	
	Sx	100	228-101	
GenEx™ Plant plus!	Mx	50	228-250	solution
	Lx	20	228-320	

### GeneAll<sup>®</sup> DirEx<sup>™</sup> series

for preperation of PCR-template without extraction

DirEx™	100	250-101	solution
DirEx <sup>™</sup> <i>Fast-</i> Tissue	96 T	260-011	solution
DirEx <sup>™</sup> <i>Fast</i> -Cultured cell	96 T	260-021	solution
DirEx <sup>™</sup> <i>Fast-</i> Whole blood	96 T	260-03 I	solution
DirEx <sup>™</sup> <i>Fast</i> -Blood stain	96 T	260-04 I	solution
DirEx <sup>™</sup> <i>Fast</i> -Hair	96 T	260-05 I	solution
DirEx <sup>™</sup> <i>Fast</i> -Buccal swab	96 T	260-06 I	solution
DirEx <sup>™</sup> <i>Fast</i> -Cigarette	96 T	260-07 I	solution

GeneAll <sup>®</sup>	RNA	series	for preperation	of total	RNA
----------------------	-----	--------	-----------------	----------	-----

RiboEx <sup>TM</sup>	mini	100	301-001	solution
		200	301-002	solution
Hybrid-R <sup>™</sup>	mini	100	305-101	spin
Hybrid-R <sup>™</sup> Blood RN.	A mini	50	3 5- 50	spin
Hybrid-R <sup>™</sup> miRNA	mini	50	325-150	spin
RiboEx <sup>™</sup> LS	mini	100	302-001	solution
RIDOEX LS	mini	200	302-002	SOlULION
Riboclear™	mini	50	303-150	spin
Riboclear <sup>™</sup> <i>plus!</i>	mini	50	3 3- 50	spin
Ribospin <sup>™</sup>	mini	50	304-150	spin
Ribospin <sup>™</sup> II	mini	50	3 4- 50	anin
	mini	300	314-103	spin
Ribospin <sup>™</sup> vRD	mini	50	302-150	spin
Ribospin <sup>™</sup> vRD <i>plus!</i>	mini	50	3 2- 50	spin
Ribospin <sup>™</sup> vRD II	mini	50	322-150	spin
Ribospin <sup>™</sup> Plant	mini	50	307-150	spin
Ribospin <sup>™</sup> Seed / Fruit	mini	50	317-150	spin
Allspin <sup>™</sup>	mini	50	306-150	spin
RiboSaver™	mini	100	351-001	solution

Products	Scale	Size	Cat. No.	Туре
GeneAll® AmpC	<b>NE<sup>™</sup></b> for	r PCR ar	nplification	
		250 L	501-025	
Taq DNA polymera	se	500 L	501-050	(2.5 U/µℓ)
		1,000 L	501-100	
		250 L	502-025	
lpha -Taq DNA polym	erase	500 L	502-050	(2.5 U/ <b>µℓ</b> )
		1,000 L	502-100	
		250 L	504-025	
lpha-Pfu DNA polyme	erase	500 L	504-050	(2.5 U/µℓ)
		1,000 L	504-100	
		250 L	505-025	
Fast-Pfu DNA polymerase		500 L	505-050	(2.5 U/µℓ)
polymerase		I,000 L	505-100	
		250 L	531-025	
Hotstart Taq DNA polymerase		500 L	531-050	(2.5 U/µℓ)
polymerase		I,000 L	531-100	
		20 µl	521-200	h an a la Roma al
Teo Decesio	0( + + + -	50 µl	521-500	- lyophilized
Taq Premix	96 tubes	20 µl	526-200	solution
		50 µl	526-500	solution
		20 µl	522-200	1 1 22 1
or Tra Decembr	0(++++-	50 µl	522-500	- lyophilized
lpha -Taq Premix	96 tube:	20 µl	527-200	and stress
		50 µl	527-500	- solution
		20 µl	525-200	1.2
HS-Taq Premix	96 tubes	s 50 μl	525-500	solution
		20 µl	520-200	lyophilized
lpha -Pfu Premix	96 tubes	s 50 μl	523-500	solution
Taq Premix (w/o dye)	96 tubes	s 20 μl	524-200	lyophilized
dNTPs mix		500 µl	509-020	2.5 mM ea
dNTPs set (set of dATP, dCTP, dGTP ar	nd dTTP)	l ml x 4 tubes	509-040	100 mM

Products	Scale	Size	Cat. No.	Туре
----------	-------	------	----------	------

### GeneAll<sup>®</sup> AmpMaster<sup>™</sup> for PCR amplification

Tra Mastaria	0.5 ml x 2 tubes	541-010	solution
Taq Master mix	0.5 ml x 10 tubes	541-050	solution
	0.5 ml x 2 tubes	542-010	solution
lpha -Taq Master mix	0.5 ml x 10 tubes	542-050	solution
	0.5 ml x 2 tubes	545-010	solution
HS-Taq Master mix	0.5 ml x 10 tubes	545-050	solution
lpha -Pfu Master mix	0.5 ml x 2 tubes	543-010	solution
	0.5 ml x 10 tubes	543-050	solution

### GeneAll<sup>®</sup> HyperScript<sup>™</sup> for Reverse Transcription

Reverse Transcript	ase 10,000 U	601-100	solution
RT Master mix	$0.5 \ {\rm ml}  imes 2 \ {\rm tubes}$	601-710	solution
RT Master mix with oligo (dT) <sub>20</sub>	$0.5 \text{ ml}  imes 2  ext{ tubes}$	601-730	solution
RT Master mix with random hexamer	0.5 ml $ imes$ 2 tubes	601-740	solution
RT Premix	96 tubes, 20 <i>µl</i>	601-602	solution
RT Premix with oligo (dT) <sub>20</sub>	96 tubes, 20 µ <b>l</b>	601-632	solution
RT Premix with random hexamer	96 tubes, 20 µl	601-642	solution
One-step RT-PCR Master mix	$0.5 \ \mathrm{ml}  imes 2 \ \mathrm{tubes}$	602-110	solution
One-step RT-PCR Premix	96 tubes, 20 μ <b>ℓ</b>	602-102	solution
First strand Synthesis Kit	50 reaction	605-005	solution
ZymAll <sup>™</sup> RNase Inhibitor	10,000 U	605-010	solution
ZymAll <sup>™</sup> RNase Inhibitor	4,000 U	605-004	solution

### GeneAll<sup>®</sup> RealAmp<sup>™</sup> for qPCR amplification

SYBR qPCR Master	200 rxn	20 <i>µl</i>	801-020	solution
mix (2X, Low ROX)	500 rxn	20 <i>µl</i>	801-050	SOIULION
SYBR qPCR Master	200 rxn	20 <i>µl</i>	801-021	solution
mix (2X, High ROX)	500 rxn	20 µl	801-051	solution

Products	Size	Cat. No.	Туре

### GeneAll<sup>®</sup> Protein series

ProtinEx <sup>™</sup> Animal cell / tissue	100 ml	701-001	solution
PAGESTA <sup>™</sup> Reducing 5X SDS-PAGE In Sample Buffer	nl × 10 tubes	751-001	solution

### GeneAll <sup>®</sup> STEAD $i^{^{\mathrm{TM}}}$ for automatic nucleic acid puritication

12 Instrument		GST012	system
24 Instrument		GST024	system
Genomic DNA Cell / Tissue	96	401-104	kit
Genomic DNA Blood	96	402-105	kit
Total RNA	96	404-304	kit
Viral DNA / RNA	96	405-322	kit
CFC Seed DNA / RNA	96	406-C02	kit
Genomic DNA Plant	96	407-117	kit
Soil DNA	96	408-114	kit

## Note

## Note

# Note



www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro, Songpa-gu, Seoul, South Korea 05729 E-mail : sales@geneall.com

Tel : 82-2-407-0096 Fax : 82-2-407-0779

©2017 GeneAll Biotechnology, All rights reserved