Ver I.I

Ribospin[™] Seed/Fruit

TOTAL RNA PURIFICATION HANDBOOK



Customer & Technical Support

Should you have any further questions, do not hesitate to contact us. We appreciate your comments and advice.

Contact Information

www.geneall.com Tel : 82-2-407-0096 Fax : 82-2-407-0779 E-mail(Order/Sales) : sales@geneall.com E-mail(Tech. Info.) : tech@geneall.com

Visit GeneAll[®] Community

www.geneall.com www.geneall.co.kr

Trademarks

Hybrid-Q[™], Exprep[™], Exfection[™], Expin[™], Exgene[™], *GenEx[™]*, DirEx[™], Hybrid-R[™], RiboEx[™], Ribospin[™], Riboclear[™], Allspin[™], Allspin[™], RiboSaver[™], EzClear[™], EzSep[™], EzPure[™], EzPass[™], AmpONE[™], AmpMaster[™], RealAmp[™], HyperScript[™], ProtinEx[™], PAGESTA[™], STEADi[™] are trademarks of GeneAll Biotechnology Co., Ltd.


© 2017 GeneAll Biotechnology, all rights reserved.

This protocol handbook is included in :

GeneAll[®] Ribospin[™] Seed/Fruit (317-150)

Visit www.geneall.com or www.geneall.co.kr for FAQ, Q&A and more information.

Brief Protocol





INDEX

Index	05
Kit Contents	06
Materials Not Provided	
Product Specifications	
Quality Control	07
Storage Conditions	
Safety Information	
Product Disclaimer	
Prevention of RNase Contamination	08
Preparation of DNase I Solution	
Protocol Selecting Guide for Starting sample	09
Product Description	10
Protocol I. For seed and fruit	11
Protocol II. For starch-enriched grain and rhizome	14
Troubleshooting Guide	17
Appendix I.	19
Appendix II.	20
Ordering Information	21

Kit Contents

Cat. No.	317-150	Storage
Components	Quantity	Storage
No. of preparation	50	
Buffer SL	30 ml	
Buffer ML	30 ml	
Buffer RBW (concentrate) *	27 ml	
Buffer RNW (concentrate) $*$ [†]	6 ml	
Buffer DRB	5 ml	Room temperature
DNase I (lyophilized) **	240 Kunitz units	(I5~25°C)
Nuclease-free water	l 5 ml	
Mini column type F (with collection tube)	50	
$EzPure^{TM}$ filter (with collection tube)	50	
1.5 ml microcentrifuge tube	50	
Protocol Handbook	I	

* Before first use, add absolute ethanol (ACS grade or better) into Buffer RBW and RNW as indicated on the bottle

[†] Contains sodium azide as a preservative

** Refer to instruction of DNase I on page 8

Materials Not Provided

- Reagent : Absolute ethanol, β -mercaptoethanol (ACS grade or better)
- Disposable material : RNase-free pipette tips, Disposable gloves
- Equipment : Microcentrifuge, Vortex mixer, Equipment for disrupting sample

Product Specifications

Ribospin™ Seed/Fruit	
Туре	Spin
Maximum amount of starting samples	100 mg/prep
No. of preparation	50
Preparation time	~ 30 min
Maximum loading volume of mini column	750 <i>µ</i> I
Minimum elution volume	30 <i>µ</i> I

Quality Control

All components in GeneAll[®] Ribospin[™] Seed/Fruit are manufactured in strictly clean conditions, and its degree of cleanness is monitored periodically. Quality control is carried out thoroughly from lot to lot, and only the qualified kits are approved to be delivered.

Storage Conditions

All components of GeneAll[®] RibospinTM Seed/Fruit should be stored at room temperature ($15\sim25^{\circ}$ C). After reconstitution of DNase I with Nuclease-free water, the DNase I solution should be stored at -20°C in aliquots for conservation of activity or used immediately for experiments.

During shipment or storage under cool ambient condition, a precipitate can form in Buffer ML. In such a case, heat the bottle to 56°C to dissolve completely. GeneAll[®] Ribospin[™] Seed/Fruit is guaranteed until the expiration date printed on the product box.

Safety Information The buffers included in the GeneAll[®] Ribospin[™] Seed/Fruit contain irritants which are harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken when handling such materials. Always wear gloves and eye protection, and follow standard safety precautions.

Buffer ML contains chaotropes, which can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

Product Disclaimer

GeneAll[®] Ribospin[™] Seed/Fruit is for research use only, not for use in diagnostic procedure.

Prevention

of RNase Contamination RNase can be introduced accidentally during RNA purification. Wear disposable gloves always, because skin often contains bacteria and molds that can be a source of RNase contamination. Use sterile, disposable plastic wares and automatic pipettes to prevent cross-contamination of RNase from shared equipment.

Preparation

of DNase I Solution The DNase I is provided in a lyophilized format. It should be reconstituted thoroughly with Nuclease-free water (provided for RNA elution) before experiment.

To obtain DNase I solution, add 120 μ I Nuclease-free water to the tube containing lyophilized DNase I (240 Kuniz units), and mix carefully and gently to avoid foaming. Dissolve the DNase I thoroughly, divide it into conveniently sized aliquots, and store at -20°C. For one preparation, 2 μ I DNase I solution is required.

Protocol Selecting Guide for Starting sample

	The list of sample applied with Protocol I	The list of sample applied with Protocol II
Seeds	Capsella bursapastoris (Shepherd's purse) Ulmus davidiana var. japonica (Elm) Daucus carota (Carrot) Raphanus sativus var. sativus (Radish) Zinnia violacea (Garden zinnia) Prunus armeniaca (Apricot tree) Apium graveolens (Celery) Pastinaca sativa (Parsley) Vitis vinifera (Grape tree) Cucurbita spp. (Pumpkin) etc.	Phaseolus vulgaris (Kidney bean) Phaseolus radiatus (Mung beans) Triticum aestivum (Wheat) Zea mays (Corn) Setaria italica (Millet) etc.
Fruits	Fragaria ananassa (Strawberry) Malus domestica (Apple) Solanum lycopersicum (Tomato) Musa sapientum L. (Banana) Mangifera indica (Mango) Pyrus serotina (Pear) Citrus unshiu (Mandarin) etc.	
Rhizomes		lpomoea batatas (Sweet potato) Solanum tuberosum (Potato) Dioscorea opposita (Yam) etc.

9

Product Description

Ribospin[™] Seed/Fruit is designed for easy and convenient isolation of total RNA from difficult plant tissues such as seeds, fruits, and rhizomes. Especially, this kit can remove effectively large quantities of secondary metabolites including polysaccharides and polyphenolic compounds which can lead to inhibition of downstream application.

Ribospin[™] Seed/Fruit provides two different procedures that are available for application of various plant tissues as follows : Protocol I for seed and fruit, Protocol II for starchenriched grain and rhizome. For efficient RNA purification, this kit offers optimized lysis system according to the sample type and adopts EzPure[™] filter column to eliminate impurities simply from lysate. Moreover, contamination of genomic DNA that causes interference in RNA analysis can be excluded by on-column DNase I treatment in these procedures.

The purified RNA is suitable for use in various downstream procedures including cDNA synthesis, RT-PCR, or Northern blotting etc.

For seed and fruit

Before starting

Thaw DNase I enzyme for use on ice. Prepare DNase I reaction mixture just before step 9. (DNase I reaction mixture : Mix 2 μ I DNase I solution with 70 μ I Buffer DRB)

1. Grind sample to a fine powder completely using a mortar and pestle under liquid nitrogen. Place up to 100 mg of ground sample into a 1.5 ml microcentrifuge tube (not provided).

Quick and complete pulverization with liquid nitrogen is essential for good result in preparation. The commonly used technique for disruption is grinding with a mortar and pestle, however other method such as bead-beater or blender can be a good alternative.

2. Add 500 μ l Buffer SL, 500 μ l Buffer ML, and 10 μ l β -mercaptoethanol to the sample and vortex vigorously for 15 sec.

Buffer ML tends to congeal with starch-enriched samples such as grain and rhizome. It is hard to separate supernatant containing RNA from debris. Therefore, if the lysate solidifies after addition of Buffer ML, use "Protocol II" instead that is special procedure for RNA extraction from starch-enriched sample.

- 3. Incubate the mixture for 3 min at room temperature.
- Centrifuge the lysate at 13,000 rpm (≥10,000 x g) for 1 min and transfer 600 µl of the supernatant to an EzPure[™] filter (yellow).
- 5. Centrifuge at 13,000 rpm (\geq 10,000 x g) for 1 min and transfer 500 μ l of the pass-through to a new 1.5 ml microcentrifuge tube (not provided).

Through this step, large cell debris and most of genomic DNAs are filtered on the $EzPure^{TM}$ filter and small pellet as debris will be formed at the bottom of the collection tube. Be careful not to disturb the pellet when transferring supernatant.

6. Add 250 μ I absolute ethanol to the supernatant and mix it well by inversion.

Do not centrifuge at this step.

After addition of absolute ethanol, precipitates may be visible in the mixture which not affect RNA purification.

7. Apply all of the mixture into a mini column type F (blue ring) and centrifuge at 13,000 rpm (\geq 10,000 x g) for I min.

Transfer all solution including any precipitates on the mini column. After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

8. Add 500 μ I Buffer RBW to the mini column and centrifuge at 13,000 rpm (\geq 10,000 x g) for 30 sec.

After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

9. Apply 70 μ I DNase I reaction mixture onto the center of the mini column for gDNA digestion. Incubate for 10 min at room temperature.

To make DNase I reaction mixture, prepare 2 μ I DNase I solution with 70 μ I Buffer DRB per on extraction. DNase I is sensitive to physical damage. Therefore, do not mix vigorously. If you want to DNase I treatment in RNA eluate, skip step 9~10 and refer to "Appendix I".

|0. Add 500 μ l Buffer RBW to the mini column and centrifuge at 13,000 rpm (\geq 10,000 x g) for 30 sec.

After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

|]. Add 500 μ l Buffer RNW to the mini column and centrifuge at 13,000 rpm (\geq 10,000 x g) for 30 sec.

After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

12. Centrifuge at maximum speed for an additional 1 min to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reaction. Make sure that the membrane of column has to be dried completely.

Protocol |

13. Add 50 µl Nuclease-free water to the center of the membrane in mini column and centrifuge at 13,000 rpm (≥10,000 x g) for 1 min.

To increase the RNA concentration, reduce the volume of elution to 30 μ l. The purified RNA should be put on ice immediately for accurate analysis or stored at -70°C for long-term storage.

For starch-enriched grain and rhizome

Before starting

Thaw DNase I enzyme for use on ice. Prepare DNase I reaction mixture just before step 10. (DNase I reaction mixture : Mix 2 μ I DNase I solution with 70 μ I Buffer DRB)

[. Grind sample to a fine powder completely using a mortar and pestle under liquid nitrogen. Place up to 100 mg of ground sample into a 1.5 ml microcentrifuge tube (not provided).

Quick and complete pulverization with liquid nitrogen is essential for good result in preparation. The commonly used technique for disruption is grinding with a mortar and pestle, however other method such as bead-beater or blender can be a good alternative.

- 2. Add 500 μ I Buffer SL and 5 μ I β -mercaptoethanol to the sample and vortex vigorously for 15 sec.
- 3. Incubate the mixture for 3 min at room temperature.
- 4. Centrifuge the lysate at 13,000 rpm (\geq 10,000 x g) for 1 min and transfer 300 μ l of the supernatant to a new 1.5 ml microcentrifuge tube (not provided).
- Add 300 µl Buffer ML to the supernatant and vortex vigorously for 15 sec and transfer all of the mixture to an EzPure[™] filter (yellow).
- 6. Centrifuge at 13,000 rpm (≥10,000 x g) for 1 min and transfer 500 µl of the pass-through to a new 1.5 ml microcentrifuge tube (not provided). Through this step, large cell debris and most of genomic DNAs are filtered on the EzPure[™] filter and small pellet as debris will be formed at the bottom of the collection tube. Be careful not to disturb the pellet when transferring supernatant.
- 7. Add 250 μ l absolute ethanol to the supernatant and mix it well by inversion.

Do not centrifuge at this step.

After addition of absolute ethanol, precipitates may be visible in the mixture which not affect RNA purification.

8. Apply all of the mixture into a mini column type F (blue ring) and centrifuge at 13,000 rpm ($\geq 10,000 \times g$) for I min.

Transfer all solution including any precipitates on the mini column. After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

9. Add 500 μ l Buffer RBW to the mini column and centrifuge at 13,000 rpm (\geq 10,000 x g) for 30 sec.

After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

10. Apply 70 μ I DNase I reaction mixture onto the center of the mini column for gDNA digestion. Incubate for 10 min at room temperature.

To make DNase I reaction mixture, prepare 2 μ I DNase I solution with 70 μ I Buffer DRB per one extraction. DNase I is sensitive to physical damage. Therefore, do not mix vigorously. If you want to DNase I treatment in RNA eluate, skip step 10~11 and refer to "Appendix I".

|]. Add 500 μ l Buffer RBW to the mini column and centrifuge at 13,000 rpm (\geq 10,000 x g) for 30 sec.

After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

12. Add 500 μ I Buffer RNW to the mini column and centrifuge at 13,000 rpm (\geq 10,000 x g) for 30 sec.

After centrifugation, discard the pass-through and re-insert the mini column back into the same collection tube.

13. Centrifuge at maximum speed for an additional 1 min to remove residual wash buffer. Transfer the mini column to a new 1.5 ml microcentrifuge tube (provided).

Residual ethanol may interfere with downstream reaction. Make sure that the membrane of column has to be dried completely.

14. Add 50 µl Nuclease-free water to the center of the membrane in mini column and centrifuge at 13,000 rpm (≥10,000 x g) for 1 min.

To increase the RNA concentration, reduce the volume of elution to 30 μ l. The purified RNA should be put on ice immediately for accurate analysis or stored at -70°C for long-term storage.

Troubleshooting Guide

Facts	Possible Causes	Suggestions
Low or no yield	Incorrect use of lysis buffer	According to the sample type, the process for lysis is different. Starch-enriched samples such as grains and rhizomes have to be processed by "Protocol II" for effective RNA extraction.
	Too much starting sample	Using too much sample leads to inefficient lysis followed by poor RNA yield. Reduce the amount of starting material.
	Insufficient pulverization	For best result, sample should be disrupted completely using proper method.
	Too low RNA mass in sample	Some samples have low RNA contents. To increase the RNA concentration in eluate, reduce the volume of elution to 30 μ l.
RNA degradation	Incorrect treatment of β-mercaptoethanol during lysis	Ensure that the correct volume of β -mercap- toethanol is used in lysis buffer for RNase elimination. The effective amount of β -mercapto- ethanol is 1% of the lysis volume.
	Improper storage of extracted RNA	The purified RNA should be stored at -70°C for long-term storage. Do not store at -20°C. If possible, perform downstream application immediately for accurate analysis after RNA extraction.
	RNase contamination	To prevent RNA degradation, wear gloves during all procedure and use RNase-free products with sterile and disposable plastic ware.
	Too old starting sample	After sufficient pulverization of starting material, store the sample properly at -70°C. If possible, perform the procedure of RNA extraction immediately after disruption of sample to decrease RNA degradation.

Troubleshooting Guide

Facts	Possible Causes	Suggestions
Clogging of EzPure™ filter	Solidification of lysate	According to the sample type, the process for lysis is different. If the lysate solidifies during lysis of protocol I, the sample may contain a lot of carbohydrate and polysaccharide. For effective RNA extraction from starch-enriched samples, apply "Protocol II".
	High viscosity of lysate	RNA can be sheared in viscous lysate that causes clogging of column. Increase centrifugal g-force and time to solve clogging if necessary.
Clogging of mini column type F	Low centrifugal force	Increase g-force (<10,000 x g) and time (~3 min).
турет	Opaque or viscous binding mixture	According to the sample type, the lysate mixed with ethanol becomes opaque or viscous. It does not affect RNA purification. However, if column is clogged because of these problems, increase centrifugal g-force and time until all mixture passes through the membrane of mini column.
DNA contamination of RNA eluate	High DNA mass in sample	Some plant tissues have high DNA contents. In this case, genomic DNA can be included in RNA eluate. To reduce DNA contamination effectively, refer to the appendix I "DNase I treatment in eluate".
	Incorrect treatment of DNase I reaction mixture	For sufficient enzymatic reaction, add DNase I reaction mixture onto the center of the membrane in mini column.

DNase I treatment in eluate

Appendix I describes how to use the DNase I (included in this kit) to eliminate contaminating genomic DNA in RNA eluate. For samples containing high DNA contents, this method is strongly recommended. This procedure is more efficient than on-column DNase I treatment.

Protocol

- I. Prepare the mixture as below in a 1.5 ml microcentrifuge tube.
 - 50 µl RNA eluate
 - 5 μ l Buffer DRB
 - I μ I DNase I solution
- 2. Incubate the mixture for 10 min at room temperature.
- 3. Add | μ | 0.25 M EDTA per 50 μ | eluate.
- 4. Inactivate DNase I enzyme at 75°C for 10 min.
- * For efficient DNase I treatment and clean-up of eluated RNA, use of Riboclear[™] plus (Cat.No. 313-150) is suggested.

Related product

Product	Cat.No.	Size	Features and Benefits
Riboclear [™] Plus	313-150	50 prep	 Preparation time : ~17 min High recovery rate : ~95% Stable and consistent yield Efficient removal of genomic DNA including DNase I Concentrated RNA eluate using micro column Complete removal of salt and enzymes No use of organic solvents, no ethanol precipitation

Electrophoresis method for using formaldehyde-agarose gel (Denaturing gel method)

A denaturing agarose gel is routinely used for the assessment of the quality of extracted RNA. The RNA isolated from samples forms secondary structure via intramolecular base pairing. Therefore, it is very difficult to analyze the result of electrophoresis because of migrating inaccuracy. However, the formaldehydeagarose gel denatures the secondary structure of RNA, making accurate migration.

To confirm the RNA band after electrophoresis, the gel should be transferred to a UV transilluminator. Mainly, two RNA bands are shown. If they are intact, the RNA bands should be sharp and the intensity of upper band should be about twice compared to that of the lower band.

Preparation of denaturing gel

- I. Put I g agarose in 72 ml water and heat to dissolve thoroughly.
- 2. Cool to 60°C.
- 3. Add 10 ml 10 X MOPS buffer, 18 ml 37% formaldehyde, and 1 μ l 10 mg/ml ethidium bromide (EtBr).
- 4. Mix well then pour the gel into the gel tray and cool to solidify it.
- 5. Transfer the solidified gel from tray to tank, and add enough 1 X MOPS running buffer to cover the gel.

Preparation of RNA sample

- I. Make the mixture.
 - -? μ I RNA (up to 20 μ g)
 - 2μ I 10 X MOPS electrophoresis buffer
- Incubate the mixture for 15 min at 65°C.
- 3. Chill the sample for 5 min in ice.
- **4.** Add $2 \mu I 10 \times$ formaldehyde gel-loading dye to the mixture.
- 5. Load the mixture in a denaturing gel which is covered with a sufficient I X MOPS electrophoresis buffer.
- 6. Run the gel and confirm the RNA band on transilluminator. Occasionally, destaining gel in dH_2O for several hours may be needed to increase the visibility of the RNA band.

Composition of buffers

10X MOPS buffer	10X formaldehyde gel-loading dye
- 0.2 M MOPS	- 50% glycerol
- 20 mM sodium acetate	- 10 mM EDTA
- 10 mM EDTA	- 0.25% (w/v) bromophenol blue
- pH to 7.0 with NaOH	- 0.25% (w/v) xylene cyanol FF

* Caution

- 4 μ l formaldehyde

- $10 \,\mu$ l formamide

When handling of formaldehyde-agarose gel, always use gloves and eye protector to avoid contact with skin and eyes. Especially, formaldehyde and ethidium bromide (EtBr) should be handled in a fume hood.

Ordering Information

	Scale	Size	Cat. No.	Туре	Products	Scale	Size	Cat. No.	Туре
GeneAll® Hybrid	I-Q™ fo	r rapid p	reparation of	plasmid DNA	GeneAll® Exgene	e [™] for is	olation o	f total DNA	
Plasmid Rapidprep	Plasmid Rapidprep	50	100-150			mini	100	105-101	spin /
	mini	200	100-102	spin	in -	rriini	250	105-152	vacuum
					Blood SV	Midi	26	105-226	spin /
GeneAll® Expre f	b™ for p	reparatic	n of plasmid	DNA	BIOOD 3V		100	105-201	vacuum
		50	101-150	spin /		MAXI	10	105-310	spin /
mini	mini	200	101-102	vacuum		1 IAV	26	105-326	vacuum
Plasmid SV		26	101-226			mini	100	106-101	spin /
Plasmid SV	Midi	50	101-250	spin / vacuum	Cell SV		250	106-152	vacuum
		100	101-201	Vacuum	Cell SV	MAXI	10	106-310	spin /
GeneAll® Exfect	ion TM					1 IAV	26	106-326	vacuum
for prepa	aration of	transfect	ion-grade pla	smid DNA		mini	100	108-101	spin /
		50	- 50	spin /			250	108-152	vacuum
Plasmid LE	mini	200	- 02	vacuum	Clinic SV	Midi	26	108-226	spin /
(Low Endotoxin)	Mill	26	-226	spin /	Clinic SV		100	108-201	vacuum
	Midi	100	-20	vacuum		MAXI	10	108-310	spin /
Plasmid EF	Midi	20	2 -220				26	108-326	vacuum
(Endotoxin Free)	1*IIDI	100	2 -20	spin		0	50	118-050	spin
						mini	100	7- 0	spin /

GeneAll® Expin '	" for pur	ification	of fragment D	NA			250	117-152	vacuum
		ification 50	of fragment E 102-150		Plant SV		250 26	7- 52 7-226	vacuum spin /
GeneAll® <i>Expin'</i> Gel SV	<i>for pur</i> mini	-		spin / vacuum	Plant SV	Midi			
Gel SV	mini	50	102-150	spin / vacuum	Plant SV	Midi	26	7-226	spin /
•		50 200	102-150 102-102	spin /	Plant SV		26 100	7-226 7-20	spin / vacuum
Gel SV PCR SV	mini	50 200 50	102-150 102-102 103-150	spin / vacuum spin / vacuum	Plant SV	Midi	26 100 10	7-226 7-20 7-3 0	spin / vacuum spin /
Gel SV	mini	50 200 50 200	102-150 102-102 103-150 103-102	spin / vacuum spin /		Midi MAXI	26 100 10 26	7-226 7-20 7-310 7-326	spin / vacuum spin / vacuum
PCR SV CleanUp SV	mini mini mini	50 200 50 200 50	102-150 102-102 103-150 103-102 113-102 113-102	spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini	Midi MAXI mini	26 100 10 26 50	7-226 7-20 7-310 7-326 14-150	spin / vacuum spin / vacuum spin
Gel SV PCR SV	mini	50 200 50 200 50 50 200	102-150 102-102 103-150 103-102 113-150	spin / vacuum spin / vacuum spin /	Soil DNA mini Stool DNA mini Viral DNA / RNA	Midi MAXI mini mini	26 100 26 50 50	7-226 7-20 7-310 17-326 14-150 15-150	spin / vacuum spin / vacuum spin spin spin
Gel SV PCR SV CleanUp SV	mini mini mini	50 200 50 200 50 200 50	102-150 102-102 103-150 103-102 113-102 113-102 112-150	spin / vacuum spin / vacuum spin / vacuum spin /	Soil DNA mini Stool DNA mini	Midi MAXI mini mini	26 100 26 50 50 50	117-226 117-201 117-310 117-326 114-150 115-150 128-150	spin / vacuum spin / vacuum spin spin
Gel SV PCR SV CleanUp SV Combo GP	mini mini mini mini	50 200 50 200 50 200 50 200 50 200	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102	spin / vacuum spin / vacuum spin / vacuum spin /	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA	Midi MAXI mini mini mini	26 100 26 50 50 50 50 250	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-152	spin / vacuum spin / vacuum spin spin spin spin
Gel SV PCR SV CleanUp SV	mini mini mini mini e TM for is	50 200 50 200 50 200 50 200 50 200	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102	spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA	Midi MAXI mini mini mini	26 100 26 50 50 50 50 250	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-152	spin / vacuum spin / vacuum spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP	mini mini mini mini	50 200 50 200 50 200 50 200 50 200 0lation o	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA	spin / vacuum spin / vacuum spin / vacuum spin /	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA	Midi MAXI mini mini mini for isol	26 100 26 50 50 50 50 250	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-152	spin / vacuum spin / vacuum spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP GeneAll® Exgent	mini mini mini e TM for is mini	50 200 50 200 50 200 50 200 50 200 0 <i>lation o</i> 100 250	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA	Midi MAXI mini mini mini	26 100 26 50 50 50 50 250 250	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-152 total DNA with	spin / vacuum spin / vacuum spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP	mini mini mini mini e TM for is	50 200 50 200 50 200 50 200 50 200 0lation o 100	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226	spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll[®] GenEx	Midi MAXI mini mini mini for isol	26 100 26 50 50 50 250 250 250	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-152 total DNA wite 220-101	spin / vacuum spin / vacuum spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP GeneAll® Exgent	mini mini mini mini e TM for is mini Midi	50 200 50 200 50 200 50 200 50 200 0 0 0 0 0 0 0 0 0 0 0 0	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll[®] GenEx	Midi MAXI mini mini mini for isol Sx Lx	26 100 26 50 50 50 250 250 250 250	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-152 total DNA wite 220-101 220-105	spin / vacuum spin / vacuum spin spin spin spin spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP GeneAll® Exgent	mini mini mini e TM for is mini	50 200 50 200 50 200 50 200 50 200 0 0 0 0 0 0 0 0 0 0 0 0	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll[®] GenEx	Midi MAXI mini mini mini for isol	26 100 26 50 50 50 250 250 250 250 250 250 100 500 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-152 total DNA wite 220-101 220-105 220-301	spin / vacuum spin / vacuum spin spin spin spin thout spin co
Gel SV PCR SV CleanUp SV Combo GP GeneAll® Exgent	mini mini mini mini e TM for is mini Midi	50 200 50 200 50 200 50 200 50 200 20	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll[®] GenEx GenEx TM Blood	Midi MAXI mini mini mini for isol Sx Lx	26 100 26 50 50 50 250 250 250 100 100 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-150 138-152 total DNA wite 220-101 220-105 220-301 221-101	spin / vacuum spin / vacuum spin spin spin spin spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP GeneAll® Exgent	mini mini mini mini e TM for is mini Midi	50 200 50 200 50 200 50 200 50 200 0 0 0 0 0 0 0 0 0 0 0 0	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll[®] GenEx GenEx TM Blood	Midi MAXI mini mini mini for isol Sx Lx Sx Lx	26 100 10 26 50 50 50 250 250 250 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-150 138-152 total DNA wite 220-101 220-105 220-301 221-101 221-105	spin / vacuum spin / vacuum spin spin spin spin spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP CeneAll® Exgent Tissue SV	mini mini mini mini mini Midi MAXI mini	50 200 50 200 50 200 50 200 0 0 0 0 0 0 0 0 0 0 0 0	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll[®] GenEx GenEx TM Blood	Midi MAXI mini mini mini for isol Sx Lx Sx	26 100 26 50 50 50 250 250 250 250 100 100 100 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-152 total DNA wite 220-101 220-105 220-301 221-101 221-105 221-301	spin / vacuum spin / vacuum spin spin spin spin spin spin solution solution
Gel SV PCR SV CleanUp SV Combo GP GeneAll® Exgent	mini mini mini mini mini mini Midi MAXI	50 200 50 200 50 200 50 200 50 200 0 0 0 0 0 0 0 0 0 0 0 0	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 ftotal DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152 109-226	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll® GenEx GenEx TM Blood GenEx TM Cell	Midi MAXI mini mini mini for isol Sx Lx Sx Lx	26 100 26 50 50 50 250 250 250 100 100 100 100 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-152 total DNA wite 220-101 220-105 220-301 221-101 221-105 221-301 222-101	spin / vacuum spin / vacuum spin spin spin spin spin spin spin spin
Gel SV PCR SV CleanUp SV Combo GP CeneAll® Exgent Tissue SV	mini mini mini mini mini Midi MAXI mini	50 200 50 200 50 200 50 200 50 200 0 0 0 0 0 0 0 0 0 0 0 0	102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA FFPE Tissue DNA GeneAll® GenEx GenEx TM Blood GenEx TM Cell	Midi MAXI mini mini mini for isol Sx Lx Sx Lx Sx	26 100 26 50 50 50 250 250 250 100 100 100 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 138-150 138-152 total DNA wite 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin / vacuum spin / vacuum spin spin spin spin spin spin spin solution solution solution solution

Products	Scale	Size	Cat. No.

GeneAll[®] GenExTM for isolation of total DNA

	Sx	100	227-101	
GenEx [™] Plant	Mx	100	227-201	solution
	Lx	100	227-301	
	Sx	100	228-101	
GenEx [™] Plant plus!	Mx	50	228-250	solution
	Lx	20	228-320	

Туре

GeneAll[®] *DirEx™* series

for preperation of PCR-template without extraction						
DirEx [™]	100	250-101	solution			
DirEx [™] <i>Fast-</i> Tissue	96 T	260-011	solution			
DirEx [™] <i>Fast</i> -Cultured cell	96 T	260-021	solution			
DirEx [™] <i>Fast-</i> Whole blood	96 T	260-03 I	solution			
DirEx [™] <i>Fast</i> -Blood stain	96 T	260-041	solution			
DirEx [™] <i>Fast-</i> Hair	96 T	260-051	solution			
DirEx [™] <i>Fast-</i> Buccal swab	96 T	260-061	solution			
DirEx [™] <i>Fast</i> -Cigarette	96 T	260-071	solution			

GeneAll[®] RNA series for preperation of total RNA

defication and the		1 1	1	
RiboEx™	mini	100	301-001	solution
NUUEX	mini	200	301-002	SOlULION
Hybrid-R [™]	mini	100	305-101	spin
Hybrid-R [™] Blood RN/	۹ mini	50	315-150	spin
Hybrid-R [™] miRNA	mini	50	325-150	spin
RiboEx [™] LS	mini	100	302-001	solution
NIDOEX L3	TT IIT II	200	302-002	SOIULION
Riboclear™	mini	50	303-150	spin
Riboclear [™] plus!	mini	50	3 3- 50	spin
Ribospin™	mini	50	304-150	spin
Ribospin™II	mini	50	3 4- 50	spin
Ribospin II	TT IIT II	300	3 4- 03	
Ribospin [™] vRD	mini	50	302-150	spin
Ribospin [™] vRD <i>plus!</i>	mini	50	3 2- 50	spin
Ribospin [™] vRD II	mini	50	322-150	spin
Ribospin [™] Plant	mini	50	307-150	spin
Ribospin [™] Seed / Fruit	mini	50	317-150	spin
Allspin™	mini	50	306-150	spin
RiboSaver™	mini	100	351-001	solution

Products	Scale	Size	Cat.	No.	Туре
GeneAll® AmpO	NE[™] for	r PCR a	mplific	ation	
		250 L		I-025	
Taq DNA polymera	se	500 L	J 50	I-050	(2.5 U/µℓ)
		I,000 L	J 50	- 00	
		250 L	J 50	2-025	
lpha-Taq DNA polyme	erase	500 L	J 50	2-050	(2.5 ∪/ µl)
		I,000 L	J 50	2-100	
		250 L	J 50	4-025	
lpha -Pfu DNA polyme	erase	500 L	J 50	4-050	(2.5 U/ µℓ)
		I,000 L	J 50	4-100	
		250 L	J 50	5-025	
Fast-Pfu DNA polymerase		500 L	J 50	5-050	(2.5 U/µℓ)
polymerase		I,000 U	J 50	5-100	
		250 L	J 53	I-025	
Hotstart Taq DNA polymerase		500 L	J 53	1-050	(2.5 U/ µℓ)
polymenase		I ,000 l	J 53	- 00	
	96 tubes	20 µl	52	I-200	han hiller d
T D '		50 μ l	52	I-500	- lyophilized
Taq Premix		20 μl	52	6-200	1.0
		50 µl	52	6-500	- solution
		20 µl	52	2-200	1 1 22 1
ar T D '	0()	50 μ l	52	2-500	- lyophilized
lpha -Taq Premix	96 tube:	20 μl	52	7-200	and streng
		50 μ l	52	7-500	- solution
		20 µl	52	5-200	1.2
HS-Taq Premix	96 tube	s 50 μl	52	5-500	- solution
		20 µl	52	0-200	lyophilized
lpha -Pfu Premix	96 tube	s 50 μθ	52	3-500	solution
Taq Premix (w/o dye)	96 tube:	s 20 μl	52	4-200	lyophilized
dNTPs mix		500 µl	50	9-020	2.5 mM eac
dNTPs set (set of dATP, dCTP, dGTP ar	nd dTTP)	l ml x 4 tubes	50	9-040	100 mM

Products	Scale	Size	Cat. No.	Туре

Products	Size	Cat. No.	Туре

GeneAll[®] AmpMaster[™] for PCR amplification

Tag Mastar pair	0.5 ml x 2 tubes	541-010	solution
Taq Master mix	0.5 ml x 10 tubes	541-050	solution
C Tag Mastar min	0.5 ml x 2 tubes	542-010	solution
lpha -Taq Master mix	0.5 ml x 10 tubes	542-050	solution
	0.5 ml x 2 tubes	545-010	solution
HS-Taq Master mix	0.5 ml x 10 tubes		solution
a Dí Maini	0.5 ml x 2 tubes	543-010	solution
lpha -Pfu Master mix	0.5 ml x 10 tubes	543-050	solution

GeneAll[®] HyperScript[™] for Reverse Transcription

Reverse Transcript	ase 10,000 U	601-100	solution
RT Master mix	$0.5 \ {\rm ml} imes 2 \ {\rm tubes}$	601-710	solution
RT Master mix with oligo (dT) ₂₀	$0.5 \text{ ml} imes 2 ext{ tubes}$	601-730	solution
RT Master mix with random hexamer	0.5 ml $ imes$ 2 tubes	601-740	solution
RT Premix	96 tubes, 20 µl	601-602	solution
RT Premix with oligo (dT) ₂₀	96 tubes, 20 μ ℓ	601-632	solution
RT Premix with random hexamer	96 tubes, 20 μ ℓ	601-642	solution
One-step RT-PCR Master mix	0.5 ml $ imes$ 2 tubes	602-110	solution
One-step RT-PCR Premix	96 tubes, 20 μ ℓ	602-102	solution
First strand Synthesis Kit	50 reaction	605-005	solution
ZymAlI [™] RNase Inhibitor	10,000 U	605-010	solution
ZymAll [™] RNase Inhibitor	4,000 U	605-004	solution

GeneAll[®] RealAmp[™] for qPCR amplification

SYBR qPCR Master	200 rxn	20 µl	801-020	solution
mix (2X, Low ROX)	500 rxn	20 <i>µl</i>	801-050	SOIULION
SYBR qPCR Master	200 rxn	20 <i>µl</i>	801-021	
mix (2X, High ROX)	500 rxn	20 <i>µl</i>	801-051	solution

GeneAll[®] Protein series

ProtinEx [™] Animal cell / tissue	100 ml	701-001	solution
PAGESTA TM Reducing 5X SDS-PAGE I mI × Sample Buffer	10 tubes	751-001	solution

GeneAll® STEAD $i^{ imes}$ for automatic nucleic acid puritication

GST012	system
GST024	system
401-104	kit
402-105	kit
403-106	kit
404-304	kit
405-322	kit
406-C02	kit
	GST024 401-104 402-105 403-106 404-304 405-322



www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro, Songpa-gu, Seoul, South Korea 05729 E-mail : sales@geneall.com

Tel : 82-2-407-0096 Fax : 82-2-407-0779

©2017 GeneAll Biotechnology, All rights reserved