

# Comparison of the Pathogen Purification Performance of Exgene™ Viral DNA/RNA to Other Commercial Kits

## Experimental Conditions

### Materials Required

- Exgene™ Viral DNA/RNA (128-150)
- 5 ml conical screw tube (for sample preparation)
- Vortex mixer
- Microcentrifuge ( $\leq 15,000 \times g$ )
- Pipette & sterile pipette tips
- Suitable protector (e.g., lab coat, disposable gloves, goggles, etc.)
- Ice (For maintenance normal state and freeze-thaw of Carrier RNA solution)

### Sample Information

- Sample type : Cultured virus and bacteria
  - Infectious bronchitis virus (IBV,  $10^{3.5} \text{EID}_{50}$ )
  - Rabies virus (RV,  $10^{3.5} \text{LD}_{50}$ )
  - Japanese encephalitis virus (JEV,  $10^{5.0} \text{TICD}_{50}$ )
  - *Mycoplasma gallisepticum* (MG,  $1 \times 10^5 \text{CCU}$ )
- Extraction conditions
  - Sample amount : 200  $\mu\text{l}$
  - Elution volume : 100  $\mu\text{l}$
  - Extraction protocol : Viral\_Normal (operation time : 29' 35")

### Sample Preparation

1. Mix the all cultured viruses and bacteria medium to 5 ml conical tube and extract the 200  $\mu\text{l}$  samples from the mixture.
2. One sample is according to Exgene™ Viral DNA/RNA protocol, the other samples are according to manual method of viral DNA/RNA extraction kits each from two different suppliers for comparison.

## Protocol

### Exgene™ Viral DNA/RNA Extraction Kit brief protocol

\* For more details and methods, please refer to [the handbook of Exgene™ Viral DNA/RNA](#).

1. Add 10  $\mu\text{l}$  of Proteinase K solution (20 mg/ml) to 1.5 ml microcentrifuge tube.
2. Transfer the 200  $\mu\text{l}$  of mixed samples and add 200  $\mu\text{l}$  of Buffer BL to the tube.
3. Add 7  $\mu\text{l}$  of Carrier RNA solution (1  $\mu\text{g}/\mu\text{l}$ ) to the tube and mix thoroughly by vortexing for 10 sec.
4. Incubate the tube at 56°C for 10 min and spin down briefly to remove any drops.
5. Add 400  $\mu\text{l}$  of Buffer RB1 to the tube and vortex for 10 sec.
6. Transfer the mixture to a Column Type S and centrifuge at  $\geq 10,000 \times g$  for 1 min at room temperature. Discard the pass-through and reinsert.
7. Add 500  $\mu\text{l}$  of Buffer BW to the column and centrifuge at  $\geq 10,000 \times g$  for 1 min at room temperature. Discard the pass-through and reinsert.
8. Add 700  $\mu\text{l}$  of Buffer TW to the column and centrifuge at  $\geq 10,000 \times g$  for 1 min at room temperature. Discard the pass-through and reinsert.
9. Centrifuge at full speed for 1 min at room temperature and transfer the column to a new 1.5 ml microcentrifuge tube.
10. Add 20~50  $\mu\text{l}$  of Nuclease-free water to the center of the membrane in the column. Stand for 1 min and centrifuge at  $\geq 10,000 \times g$  for 1 min at room temperature.

## Result





Figure 1. Comparison of C<sub>t</sub> value of DNA template extracted from the *Mycoplasma gallisepticum* (MG).

The DNA templates were extracted from the *Mycoplasma gallisepticum* (MG) using Exgene™ Viral DNA/RNA (red line) and viral DNA/RNA extraction kits (manual method) each from two different suppliers (blue and green line). Eluted DNA templates were analyzed with a TaqMan-based real-time PCR assay using CFX-96.

- Red line : GeneAll® Exgene™ Viral DNA/RNA
- Blue line : Supplier A viral DNA/RNA extraction kit
- Green line : Supplier B viral DNA/RNA extraction kit
- PCR instrument : CFX-96 (1855201)
- qPCR kit : Probe qPCR Mix (RR391A)
- Target gene : None specific


© 2022 GeneAll Biotechnology, all rights reserved.

## Comparison of the Pathogen Purification Performance of Exgene™ Viral DNA/RNA to Other Commercial Kits




**Figure 2.** Analysis of extraction sensitivity on serial dilutions of *Mycoplasma gallisepticum* (MG). The DNA templates were extracted from a 10-fold serial dilution of *Mycoplasma gallisepticum* (MG) using Exgene™ Viral DNA/RNA. All eluates were analyzed with a TaqMan-based real-time PCR assay using CFX-96.

- N : Negative control (Nuclease-free water)
- PCR instrument : CFX-96 (1855201)
- qPCR kit : Probe qPCR Mix (RR391AT)
- Target gene : None specific



**Figure 4.** Analysis of extraction sensitivity on serial dilutions of Infectious bronchitis virus (IBV). The RNA template were extracted from a 10-fold serial dilution of Infectious bronchitis virus (IBV) using Exgene™ Viral DNA/RNA. All elutes were analyzed with TaqMan-based one-step RT-qPCR assay using CFX-96.


- N : Negative control (Nuclease-free water)
- PCR instrument : CFX-96 (1855201)
- RT-qPCR kit : 2X 1 Step RT-qPCR Master Mix [for probe] (QRT1-XV-100R)
- Target gene : None specific



**Figure 3.** Comparison of C<sub>t</sub> value of DNA template extracted from the Infectious bronchitis virus (IBV).

The RNA templates were extracted from the Infectious bronchitis virus (IBV) using Exgene™ Viral DNA/RNA (red line) and viral DNA/RNA extraction kits (manual method) each from two different suppliers (blue and green line). Eluted RNA templates were synthesized to cDNA with reverse transcription; and then analyzed with TaqMan-based one-step RT-qPCR assay using CFX-96.

- Red line : GeneAll® Exgene™ Viral DNA/RNA
- Blue line : Supplier A viral DNA/RNA extraction kit
- Green line : Supplier B viral DNA/RNA manual extraction kit
- PCR instrument : CFX-96 (1855201)
- RT-qPCR kit : 2X 1 Step RT-qPCR Master Mix [for probe] (QRT1-XV-100R)
- Target gene : None specific




**Figure 5.** Analysis of extraction sensitivity on dilutions of Rabies virus (RV).

The RNA template were extracted from a 10-fold serial dilution from 10<sup>0</sup> to 10<sup>-5</sup> of known positive Rabies virus (RV) samples using Exgene™ Viral DNA/RNA. All eluates were analyzed using conventional reverse transcription PCR (one-step RT-PCR) assay.

- M : GENESTA™ 250 bp DNA ladder (GA-025)
- N : Negative control (Nuclease-free water)
- Target gene (PCR product size) : Jecom (100 bp)
- PCR instrument : MultiGene™ Optimax thermal cycler (TC9610, Supplier : L)
- RT-PCR kit : HyperScript™ One-Step RT-PCR Master Mix, 0.5 ml x 2 (602-125)
- Electrophoresis conditions : 1.2% agarose, 110 V, 30 min, 10 µl

# Comparison of the Pathogen Purification Performance of Exgene™ Viral DNA/RNA to Other Commercial Kits



**Figure 6. Analysis of extraction sensitivity on dilutions of Japanese encephalitis virus (JEV).** The RNA template were extracted from a 10-fold serial dilution from 10<sup>0</sup> to 10<sup>-4</sup> of known positive Japanese encephalitis virus (JEV) samples using Exgene™ Viral DNA/RNA. All eluates were analyzed using conventional reverse transcription PCR (one-step RT-PCR) assay.

- M : GENESTA™ 250 bp DNA ladder (GA-025)
- N : Negative control (Nuclease-free water)
- Target gene (PCR product size) : omRABV (192 bp)
- PCR instrument : MultiGene™ Optimax thermal cycler (TC9610, Supplier : L)
- RT-PCR kit : HyperScript™ One-Step RT-PCR Master Mix, 0.5 ml x 2 (602-125)
- Electrophoresis conditions : 1.2% agarose, 110V , 30 min, 10 µl