

# High-quality and Highly Scalable Cell DNA Extraction using GenEx™ Cell, the Solution-type Genomic DNA Extraction Kit

## Experimental Conditions

### Materials Required

- GenEx™ Cell Sx (221-101)
- 200 U of lyticase or 20 U of zymolase (for yeast cell lysis)
- 30 mg/ml of lysozyme or 300 µl/ml of lysostaphin (for gram-positive cell lysis)
- Ice (for incubation or maintaining the normal state of the enzyme solution and Proteinase K solution)
- Microcentrifuge tube
- Microcentrifuge (≤15,000 x g)
- Vortex mixer
- Heating block
- 70% ethanol
- Isopropanol (≥99.5%, C<sub>3</sub>H<sub>8</sub>O, CAS No. 67-63-0)
- Pipette & sterile pipette tips
- Suitable protector (e.g., lab coat, disposable gloves, goggles, etc.)

### Sample Information

- 5 x 10<sup>6</sup> of K562 (human erythroleukemia cell line)
- 2 x 10<sup>9</sup> of DH5α (gram-negative bacteria)
- 2 x 10<sup>9</sup> of *Lactobacillus* (gram-positive bacteria)
- 5 x 10<sup>7</sup> of yeast

## Protocol

Enzymatic pre-treatment of lysozyme or lysostaphin is required for gram-positive and yeast cell DNA extraction.

Gram-positive bacteria : manual handbook of Exgene™ Cell SV mini (Step 1~2 of Protocol L, page 38) was observed.

Yeast : manual handbook of Exgene™ Cell SV mini (Step 1~4 of Protocol M, page 40) was observed.

For more details, [please refer to handbook of Exgene™ Cell SV mini and GenEx™ Cell Sx.](#)

## Sample Preparation and protocol

### K562

Manual handbook of Protocol E (page 24~25) was observed.

### DH5α

Manual handbook of Protocol G (page 28~29) was observed.

### *Lactobacillus*

Following enzymatic pre-treatment, manual handbook was observed starting with Step 3 of Protocol G (page 28).

### Total yeast

Following enzymatic pre-treatment, manual handbook was observed starting with Step 2 of Protocol G (page 28).


## Result

Kit	GenEx™ Cell Sx			CV (%)
	Yield (ng/µl)	A <sub>260/280</sub>	A <sub>260/230</sub>	
K562	15.62	1.98	2.15	1.10
	15.51	2.01	2.21	
	15.29	1.98	2.18	
DH5α	14.92	1.98	2.19	0.85
	15.11	1.99	2.20	
	14.86	2.01	1.99	
<i>Lactobacillus</i>	12.26	1.98	2.19	3.78
	12.68	2.02	2.18	
	13.21	1.98	2.20	
Yeast	6.11	2.01	2.11	2.77
	6.35	1.98	2.19	
	6.02	2.03	2.20	

**Table 1. Result of yield, purity and CV (coefficient of variation) of DNA extracted from 4 different samples using GenEx™ Cell Sx.**

The DNA were extracted from four different samples using GenEx™ Cell. All eluates were analyzed with an absorbance using NanoDrop™ 2000. The absorbance was performed in triplicated. The yield and CV values were calculated based on the measured absorbance values.

# High-quality and Highly Scalable Cell DNA Extraction using GenEx™ Cell, the Solution-type Genomic DNA Extraction Kit



**Figure 1. PCR amplification of DNA extracted from four different samples using GenEx™ Cell Sx.**

PCR reaction was performed in triplicate with extracted DNA from four different samples using GenEx™ Cell Sx. Eluted PCR products were analyzed with gel-electrophoresis using ethidium bromide staining.

- **PCR primer**

- A : K562 cells : human GAPDH primer
- B : Total yeast : Scer primer
- C : DH5α : bacteria universal primer
- D : *Lactobacillus* : uvrC primer

- **DNA ladder**

- Lane M : DNA ladder
- A, B : GENESTA™ 100 bp DNA ladder (GA-010)
- C, D : GENESTA™ 250 bp DNA ladder (GA-025)

- **PCR instrument and kit**

- MultiGene™ Optimax Thermal Cycler (TC9610, Supplier : L)
- 2X Taq PCR Master Mix (TAQ-OV-500R, Supplier : M)