Ver I.O



**2010 CELL SV MAXI** 

STOOD SA WIDI/WYXI

gyggng

**DNA PURIFICATION HANDBOOK** 



### **Customer & Technical Support**

Do not hesitate to ask us any question.

We thank you for any comment or advice.

### Contact us at

www.geneall.com

Tel: 82-2-407-0096

Fax: 82-2-407-0779

E-mail(Order/Sales): sales@geneall.com

E-mail(Tech. Info.): tech@geneall.com

### Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

### **Trademarks**

AmpONE<sup>TM</sup>, Exfection<sup>TM</sup>, Exgene<sup>TM</sup>, Expin<sup>TM</sup>, Exprep<sup>TM</sup>, EzClear<sup>TM</sup>, EzSep<sup>TM</sup>,  $GenEx^{TM}$ , Hybrid-Q<sup>TM</sup>, RiboEx<sup>TM</sup>, Ribospin<sup>TM</sup> are trademarks of GeneAll Biotechnology co., Itd.

© 2016 GeneAll Biotechnology, all right reserved.

This protocol handbook is included in:

GeneAll® Exgene™ Blood SV Midi (105-226,105-201) GeneAll® Exgene™ Blood SV MAXI (105-310, 105-326)

GeneAll® Exgene™ Clinic SV Midi (108-226, 108-201)

GeneAll® Exgene™ Clinic SV MAXI (108-310, 108-326)

GeneAll® Exgene™ Cell SV MAXI (106-310, 106-326)

Visit www.geneall.com or www.geneall.co.kr for FAQ, QnA and more information.

# **INDEX**

		Kit Contents	04
		Storage Condition	05
		Introduction	06
		General Information	07
		Centrifuge for Midi/MAXI	08
		Protocols for	
Whole Blood	$0.4 \sim 1 \text{ ml}$	Blood SV Midi Clinic SV Midi	10
	$I \sim 2 \text{ ml}$	Blood SV Midi Clinic SV Midi	14
	$3 \sim 5 \text{ ml}$	Blood SV MAXI Clinic SV MAXI Cell SV MAXI	18
	$6 \sim 10  \text{ml}$	Blood SV MAXI Clinic SV MAXI Cell SV MAXI	22
Cultured cells	$\sim 2 \times 10^7$ cells	Clinic SV Midi	26
	$\sim 1 \times 10^8$ cells	Clinic SV MAXI Cell SV MAXI	30
Animal tissue	30 ∼ 100 mg	Clinic SV Midi	34
, amriar dissact	100 ~ 250 mg		38
	100 ~ 250 Hig	Clinic SV MAXI Cell SV MAXI	30
Gram (-) bacteria	$\sim 1 \times 10^{10}$ cells	Clinic SV Midi	42
	$\sim 5 \times 10^{10}$ cells	Clinic SV MAXI Cell SV MAXI	46
Gram (+) bacteria	$\sim 5 \times 10^{10} \text{ cells}$	Cell SV MAXI	50
Yeast	~ 5 x 10 <sup>8</sup> cells	Cell SV MAXI	55
		Troubleshooting Guide	60
		Ordering information	64

# CONTENTS

		Blood SV	d SV			Clinic SV	c SV		Cell SV	SV
Cat. No.	105-226	105-226 105-201 105-310 105-326	105-310	105-326		108-201	108-226 108-201 108-310 108-326	108-326	106-310 106-326	106-326
Size	Midi	Midi	MAXI	MAXI	Midi	Midi	MAXI	MAXI	MAXI	MAXI
Preps.#	26	100	10	26	26	100	10	26	10	26
Columns	26	100	10	26	26	100	10	26	10	26
Tubes	52	200	20	52	52	200	20	52	20	52
Buffer GP (ml)	ı	•	ı	1	T	ı	ı	1	09	150
Buffer YL (ml)	ı	1	ı	1	ī	ı	ı	ı	09	150
Buffer CL (ml)	ı	1	1	1	09	200	09	200	09	200
Buffer BL (ml)	80	300	150	320	80	300	150	320	150	320
Buffer BW (ml)	06	400	06	220	06	400	06	220	06	220
Buffer TW (ml)	100	200	120	360	100	200	120	360	120	360
Buffer AE (ml)	30	120	09	120	30	120	09	120	09	120
Proteinase K (mg)	09	240	48	120	09	240	48	120	48	120
PK Storage bfr. (ml)	4	14	4	7	4	14	4	7	4	7
Protocol Handbook	_	_	_	_	_	_	_	_	_	_

### **Precautions and disclaimer**

GeneAll® Exgene<sup>TM</sup> Midi and MAXI series are for research use only, and should not be used for drug, household or other unintended uses. All due care and attention should be taken in every procedure in this handbook. Please consult the Material Safety Data Sheet (MSDS) for information regarding hazards and safe handling practices.

# **Storage Condition**

All components of GeneAll® Exgene<sup>TM</sup> Midi and MAXI series should be stored at room temperature ( $15\sim25^{\circ}$ C). After reconstitution of Proteinase K with storage buffer, it should be stored under 4°C for conservation of activity. It can be stored at 4°C for I year without significant decrease in activity. But for prolonged preservation of activity, storing under -20°C is recommended.

Under cool ambient condition, a precipitate can be formed in buffer CL and/or BL. In such a case, heat the bottle above 37°C to dissolve completely. GeneAll® Exgene<sup>TM</sup> Midi and MAXI series are guaranteed until the expiration date printed on the product label.

### **Chemical Hazard**

Buffer BL and BW contain irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling. Always wear gloves and eye protector, and follow standard safety precautions.

# GeneAll® Exgene™

# **Blood/Clinic/Cell SV Midi/MAXI**

### Introduction

GeneAll® Exgene<sup>TM</sup> Midi and MAXI series provide fast and easy methods of large scale purification of total DNA from various biological samples, such as whole blood, cultured cells, tissues, bacteria and etc.

GeneAll® Exgene<sup>TM</sup> Series utilize the advanced silica-binding technology to purify total DNA sufficiently pure for many applications. Various samples are lysed in optimized buffer containing detergents and lytic enzyme. Under high salt condition, DNA in the lysate bind to silica membrane and impurities pass through membrane into a collecting tube. The membranes are washed with a series of alcohol-containing buffer to remove any traces of proteins, cellular debris and salts. Finally pure DNA is released into a clean collection tube with deionized water or low ionic strength buffer.

# **Typical yield**

Sample	Amount	Typical Yield
Human Whole blood	2 ml	20 ∼ 60 ug
Animal tissue	100 mg	20 ~ 80 ug
Cultured cells	$2 \times 10^7$ cells	40 ~ 100 ug
Bacteria cells	$1 \times 10^{10}$ cells	30 ∼ 100 ug
Yeast	$5 \times 10^8$ cells	70 ~ 200 ug

# **General Information**


- All centrifugation on procedures MUST be performed on a tabletop centrifuge which has a swinging bucket. Using fixed-angle-rotor centrifuge will cause poor yields of DNA recovery. Refer to next page for suitable centrifuge.
- Centrifugal force should be over 4,000 xg at least to get proper result.
   Low centrifugal force will cause the incomplete washing and eluting, followed by poor result.
- Unless RNase A is treated, both DNA and RNA will be co-purified.
   RNA may inhibit some downstream enzymatic reactions, but not PCR itself. If RNA-free DNA is required, RNase should be treated before addition of Buffer BL.
- $\bullet$  The size of purified DNA is up to 50kb in length, and the majority is 20  $\sim$  30 kb approximately. This fragments will be completely denatured during PCR and it can be amplified with high efficiency.
- DNA purified by GeneAll® Exgene<sup>™</sup> Midi and MAXI procedure is free of protein and other contaminants which may inhibit PCR or other enzymatic reactions.
   Purified DNA can be applied to various downstream applications, such as PCR, Southern blotting, RAPD, AFLP, RFLP and other subsequent enzymatic reactions.
   DNA can be used immediately or safely stored in buffer AE at -20°C for later use.

# **Using Swinging-bucket Centrifuge in**

# Midi/MAXI procedures

GeneAll® Exgene™ Midi and MAXI procedures require the conventional centrifuge which has a swinging-bucket rotor and ability of 4,000 ~ 5,000 xg. Use of fixed-angle rotor will cause nonuniform contact between column membrane and solutions, followed by inconsistent result. Low g-force will lead to uncomplete removal of ethanol from column membrane and to inadequate eluting. Compatible centrifuges and rotors are listed below, but you can use any other equivalent.

Company	Centrifuge	Rotor
Beckman CoulterInc. (California, USA)	Allegra X-15R Allegra 25R	Sx4750 Sx4750A TS-5.1-500
Eppendorf AG (Hamburg, Germany)	5804/5804R 5810/5810R	A-4-44
EYELA Inc. (Tokyo, Japan)	5800 5900	RS-410 RS-410M
Hanil Science Industrial Inc. (Incheon, Korea)	Union 5KR Union 55R MF-550 HA1000-6 HA1000-3	R-WS1000-6B W-WS750-6B HSR-4S WHSR-4S
Hettich AG (Kirchlengern, Germany)	Rotina 35 Rotanta 460 Rotixa 50S	1717 1724 5624



# **GeneAll<sup>®</sup> Exgene<sup>™</sup> Protocols**

Cell SV MAXI
Clinic SV Midi/MAXI
Blood SV Midi/MAXI



# PROTOCOL for 0.4 ~ I ml of whole blood





### Before experiment

Prepare the water bath to 56°C

Prepare absolute ethanol

If a precipitate has formed in buffer BL, heat to dissolve completely

- 1. Pipet 50 ul of Proteinase K solution (20 mg/ml, provided) into the bottom of a 15 ml centrifugation tube (not provided).
- 2. Add I ml of the sample to the tube and mix well.

If the sample volume is less than 1 ml, bring the volume of sample to 1 ml with PBS.

3. Add I.2 ml of Buffer BL to the tube. Vortex the tube for I5 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely.

4. Incubate at 56°C for 20 min.

During incubation, occasional vortexing of the lysate will help accelerate lysis. Longer incubation will not affect DNA recovery.

5. Add I ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

6. Transfer all of the mixture to a Midi column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of the mixture to a Midi column, be careful not to moisten the rim of a Midi column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

7. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column.

8. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the column membrane has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the Midi column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

# 9. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow this;

→ Discard the filtrate, wipe off any spillage from the thread of the 15 ml tube, and re-insert the Midi column back into the 15 ml tube. Apply 1 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the Midi column and incubate it at room temperature for 15 min.

Insufficient centrifugal force will bring on residual ethanol in Midi column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

# 10. Place the Midi column into a new 15 ml centrifugation tube (provided). Pipet 300 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge for 5 min at 4,500 xg (5,000 rpm).

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in Midi column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the Midi column membrane for optimal elution of DNA.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

If the volume of starting sample is less than 1 ml, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 100 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problems by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water is higher than 7.0.

11. For higher concentrated yield, re-load the eluate from step 10 into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

For higher total yield, add 300 ul of fresh Buffer AE or distilled water again into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

The first and second eluates can be combined or collected seperately as necessity.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

# PROTOCOL for I ~ 2 ml of whole blood





### Before experiment

Prepare the water bath to 56°C

Prepare absolute ethanol

If a precipitate has formed in buffer BL, heat to dissolve completely

- 1. Pipet 100 ul of Proteinase K solution (20 mg/ml, provided) into the bottom of a 15 ml centrifugation tube (not provided).
- 2. Add 2 ml of the sample to the tube and mix well.

If the sample volume is less than 2 ml, bring the volume of sample to 2 ml with PBS.

3. Add 2.4 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely.

4. Incubate at 56°C for 20 min.

During incubation, occasional vortexing of the lysate will help accelerate lysis. Longer incubation will not affect DNA recovery.

5. Add 2 ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

6. Transfer 4 ml of the mixture to a Midi column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of the mixture to a Midi column, be careful not to moisten the rim of a Midi column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

7. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply the remainder of the mixture, close the cap, and centrifuge for 3 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column.

While transfer of the mixture to the Midi column, be careful not to moisten the rim of Midi column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

8. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column.

9. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the column membrane has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the Midi column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

# 10. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

→ Discard the filtrate, wipe off any spillage from the thread of the 15 ml tube, and re-insert the Midi column back into the 15 ml tube. Apply 1 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the Midi column and incubate it at room temperature for 15 min.

Insufficient centrifugal force will bring on residual ethanol in Midi column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

# I I. Place the Midi column into a new 15 ml centrifugation tube (provided). Pipet 400 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge for 5 min at 4,500 xg (5,000 rpm).

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in Midi column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the Midi column membrane for optimal elution of DNA.

Less than 400 ul of eluate will be obtained from 400 ul of elution buffer, but this has no influence on DNA yields.

If the volume of starting sample is less than 2 ml, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 100 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problems by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water is higher than 7.0.

12. For higher concentrated yield, re-load the eluate from step 11 into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

Less than 400 ul of eluate will be obtained from 400 ul of elution buffer, but this has no influence on DNA yields.

For higher total yield, add 400 ul of fresh Buffer AE or distilled water again into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

The first and second eluates can be combined or collected seperately as necessity.

Less than 400 ul of eluate will be obtained from 400 ul of elution buffer, but this has no influence on DNA yields.

# PROTOCOL for 3 ~ 5 ml of whole blood

### Before experiment







Prepare the water bath to 65°C Prepare absolute ethanol If a precipitate has formed in buffer BL, heat to dissolve completely

- 1. Pipet 200 ul of Proteinase K solution (20 mg/ml, provided) into the bottom of a 50 ml conical tube (not provided).
- 2. Add 5 ml of the sample into the tube and mix well. If the sample volume is less than 5 ml, bring the volume of sample to 5 ml with PBS.
- 3. Add 6 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely and yield a homogenous mixture.

4. Incubate at 65°C for 20 min.

During incubation, occasional vortexing of the lysate will help accelerate lysis. Longer incubation will not affect DNA recovery.

5. Add 5 ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

6. Transfer all of the mixture to a MAXI column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

Do NOT place the MAXI columns in tilted or bottom-up position even if caps are closed.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

7. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply 7 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 50 ml tube before reinserting the MAXI column

8. Discard the filtrate, re-insert the MAXI column back into the 50 ml tube. Apply 10 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the MAXI column has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the MAXI column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

# 9. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

→ Discard the filtrate, wipe off any spillage from the thread of the 50 ml tube, and re-insert the MAXI column back into the 50 ml tube. Apply 3 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the MAXI column and incubate it for 20 min at room temperature.

Insufficient centrifugal force will bring on residual ethanol in MAXI column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

10. Place the MAXI column into a new 50 ml centrifugation tube (provided). Pipet 600 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in MAXI column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the MAXI column membrane for optimal elution of DNA.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yields.

If the volume of starting sample is less than 5 ml, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 300 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA (0.5mM) included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water (at least 7.0) before elution.

11. For higher concentrated yield, re-load the eluate from step 10 into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

For higher total yield, add 600 ul of fresh buffer AE or distilled water again into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

The first and second eluates can be combined or collected seperately as necessity.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

# PROTOCOL for 6 ~ 10 ml of whole blood

\* Due to the need of additional reagents, fewer preparations can be performed. Additional reagents can be purchased separately as accessory.







#### Before experiment

Prepare the water bath to 65°C

Prepare absolute ethanol

If a precipitate has formed in buffer BL, heat to dissolve completely

- 1. Pipet 400 ul of Proteinase K solution (20 mg/ml, provided) into the bottom of a 50 ml conical tube (not provided).
- 2. Add 10 ml of the sample into the tube and mix well.

If the sample volume is less than 10 ml, bring the volume of sample to 10 ml with PBS.

3. Add 12 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely and yield a homogenous mixture.

4. Incubate at 65°C for 20 min.

During incubation, occasional vortexing of the lysate will help accelerate lysis. Longer incubation will not affect DNA recovery.

5. Add 10 ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

6. Transfer a half of the mixture to a MAXI column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

Do NOT place the MAXI columns in tilted or bottom-up position even if caps are closed.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

7. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply the remainder of the mixture, close the cap, and centrifuge for 3 min at 2,000 xg (3,000 rpm)

Wipe off any spillage from the thread of the 50 ml tube before reinserting the MAXI column

While transfer of the mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

8. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply 7 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3.000 rpm).

Wipe off any spillage from the thread of the 50 ml tube before reinserting the MAXI column.

9. Discard the filtrate, re-insert the MAXI column back into the 50 ml tube. Apply 10 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the MAXI column has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the MAXI column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

# 10. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

→ Discard the filtrate, wipe off any spillage from the thread of the 50 ml tube, and re-insert the MAXI column back into the 50 ml tube. Apply 3 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the MAXI column and incubate it for 20 min at room temperature.

Insufficient centrifugal force will bring on residual ethanol in MAXI column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

# II. Place the MAXI column into a new 50 ml centrifugation tube (provided). Pipet I ml of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in MAXI column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the MAXI column membrane for optimal elution of DNA.

Less than I ml of eluate will be obtained from I ml of elution buffer, but this has no influence on DNA yields.

If the volume of starting sample is less than 10 ml, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 300 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA (0.5mM) included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water (at least 7.0) before elution.

12. For higher concentrated yield, re-load the eluate from step 11 into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Less than I ml of eluate will be obtained from I ml of elution buffer, but this has no influence on DNA yield.

For higher total yield, add I ml of fresh buffer AE or distilled water again into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

The first and second eluates can be combined or collected seperately as necessity.

Less than I ml of eluate will be obtained from I ml of elution buffer, but this has no influence on DNA yield.

# PROTOCOL for ~ 2 x 10<sup>7</sup> of cultured cells

### Before experiment



Prepare the water bath to 56°C and 70°C

Prepare absolute ethanol

If a precipitate has formed in buffer CL or BL, heat to dissolve completely

1. Pellet cells (up to  $2 \times 10^7$  cells) to a 15 ml microcentrifuge tube by centrifugation at 2,000 xg for 5 min.

Certain cell strains, such as PC12, are not lysed well in buffer CL. For those cells, it is helpful to perform additional freeze-thaw step several times before proceeding to next step.

2. Discard the supernatant as much as possible and resuspend thoroughly cell pellet in I ml of Buffer CL.

Pelleted cells may not be resuspended easily in buffer CL. It is helpful to resuspend the cell pellet with residual media by flickering or vortexing before the addition of buffer CL.

3. Add 100 ul of Proteinase K solution (20 mg/ml, provided). Mix thoroughly by vortexing. Incubate for 20 min at 56°C.

Vortex the lysate occasionally to accelerate during incubation. Longer incubation will not affect DNA recovery.

 (Optional:) If RNA-free DNA is required, cool the mixture to room temperature, add 20 ul of RNase solution (100 mg/ml, Cat.No.117-961), vortex to mix thoroughly, and incubate for 5 min at room temperature.

Unless RNase is treated, both DNA and RNA will be co-purified. RNA can inhibit some downstream enzymatic reactions, but not PCR itself.

# 5. Add 1.2 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely.

- 6. Incubate at 70°C for 10 min.
- 7. Add I ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

8. Transfer all of the mixture to a Midi column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of the mixture to a Midi column, be careful not to moisten the rim of a Midi column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

9. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of Buffer BW and centrifuge for 3 min at 2.000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column.

10. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the column membrane has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the Midi column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

# II. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow this;

→ Discard the filtrate, wipe off any spillage from the thread of the 15 ml tube, and re-insert the Midi column back into the 15 ml tube. Apply 1 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the Midi column and incubate it at room temperature for 15 min.

Insufficient centrifugal force will bring on residual ethanol in Midi column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

12. Place the Midi column into a new 15 ml centrifugation tube (provided). Pipet 300 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge for 5 min at 4,500 xg (5,000 rpm).

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in Midi column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the Midi column membrane for optimal elution of DNA.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

If the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 100 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problems by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water is higher than 7.0.

13. For higher concentrated yield, re-load the eluate from step 12 into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

For higher total yield, add 300 ul of fresh Buffer AE or distilled water again into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

The first and second eluates can be combined or collected seperately as necessity.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

# PROTOCOL for ~ I x 10<sup>8</sup> of cultured cells

#### Before experiment





Prepare the water bath to 65°C and 70°C

Prepare absolute ethanol

If a precipitate has formed in buffer CL or BL, heat to dissolve completely

1. Pellet cells (up to 1x10<sup>8</sup> cells) to a 50 ml microcentrifuge tube by centrifugation at 2,000 xg for 5 min.

Certain cell strains, such as PC12, are not lysed well in buffer CL. For those cells, it is helpful to perform additional freeze-thaw step several times before proceeding to next step.

2. Discard the supernatant as much as possible and resuspend thoroughly cell pellet in 5 ml of Buffer CL.

Pelleted cells may not be resuspended easily in buffer CL. It is helpful to resuspend the cell pellet with residual media by flickering or vortexing before the addition of buffer CL.

3. Add 200 ul of Proteinase K solution (20 mg/ml, provided). Mix thoroughly by vortexing. Incubate for 20 min at 65°C.

Vortex the lysate occasionally to accelerate during incubation. Longer incubation will not affect DNA recovery.

4. (Optional:) If RNA-free DNA is required, cool the mixture to room temperature, add 100 ul of RNase solution (100 mg/ml, Cat. No.117-961), vortex to mix thoroughly, and incubate for 5 min at room temperature.

Unless RNase is treated, both DNA and RNA will be co-purified. RNA can inhibit some downstream enzymatic reactions, but not PCR itself.

### 5. Add 6 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely and yield a homogenous mixture.

#### 6. Incubate at 70°C for 10 min.

### 7. Add 5 ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

# 8. Transfer all of the mixture to a MAXI column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

Do NOT place the MAXI columns in tilted or bottom-up position even if caps are closed.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

# 9. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply 7 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 50 ml tube before reinserting the MAXI column.

# 10. Discard the filtrate, re-insert the MAXI column back into the 50 ml tube. Apply 10 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5.000 rpm).

If the MAXI column has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the MAXI column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

# II. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

→ Discard the filtrate, wipe off any spillage from the thread of the 50 ml tube, and re-insert the MAXI column back into the 50 ml tube. Apply 3 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the MAXI column and incubate it for 20 min at room temperature.

Insufficient centrifugal force will bring on residual ethanol in MAXI column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

12. Place the MAXI column into a new 50 ml centrifugation tube (provided). Pipet 600 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in MAXI column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the MAXI column membrane for optimal elution of DNA.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yields.

If the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 300 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA (0.5mM) included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water (at least 7.0) before elution.

13. For higher concentrated yield, re-load the eluate from step 12 into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

For higher total yield, add 600 ul of fresh buffer AE or distilled water again into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

The first and second eluates can be combined or collected seperately as necessity.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

# G PROTOCOL for 30 ~ 100 mg animal tissue

#### Before experiment



Prepare the water bath to 56°C and 70°C

Prepare absolute ethanol

If a precipitate has formed in buffer CL or BL, heat to dissolve completely

1. Homogenize  $30 \sim 100$  mg of tissue as described in step IA, IB or IC, depending on the sample type.

It is most important to weigh the sample accurately.

If the sample is spleen tissue, up to 40 mg can be processed.

Well-disrupted sample will accelerate lysis and decrease the lysis time.

- IA. For soft tissue, such as liver or brain, put up to 100 mg of the tissue into 15 ml conical tube, add 400 ul of Buffer CL, homogenize thoroughly on ice with homogenizer, add 600 ul of Buffer CL, and vortex vigorously to homogenate well.
- IB. If a homogenizer is not available or the tissue is not soft, grind the tissue to a fine powder with liquid nitrogen in a pre-cooled mortar and pestle. Put up to 100 mg of the powdered tissue into 15 ml conical tube. Add I ml of Buffer CL and vortex for 30 sec to homogenate completely.
- IC. If neither IA nor IB is available, mince the tissue with sterile sharp blade or scalpel as small as possible. Put up to 100 mg of the tissue into a 15 ml conical tube. Add I ml of Buffer CL and pulse-vortex for 30 sec.
  - \*\* Alternatively, tissue sample can be effectively disrupted using some instruments, such as rotor-stator homogenizer or a bead-beater.

    When use these, follow the manufacture s instruction manual.

2. Add 100 ul of Proteinase K solution to the tube. Mix completely by vortexing or pipetting. Incubate at 56°C until the sample is completely lysed.

It is essential to mix the components completely for efficient lysis.

Lysis time varies from 10 min to 3 hr usually depending on the type of tissue and the homogenization method. The lysate should become translucent without any particles after complete lysis. Overnight lysis does not influence the preparation. If the sample is lysed in water bath or heating block, vortex occasionally (2  $\sim$  3 times per hour) during incubation to lyse readily. Lysis in shaking water bath, shaking incubator, or agitator would be best for efficient lysis.

3. (Optional:) If RNA-free DNA is required, cool the lysate to room temperature, add 20 ul of RNase solution (100 mg/ml, Cat.No.117-961), vortex to mix thoroughly, and incubate for 3 min at room temperature.

Unless RNase is treated, both DNA and RNA will be co-purified. RNA can inhibit some downstream enzymatic reactions, but not PCR itself.

4. Add I.2 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely.

- 5. Incubate at 70°C for 10 min.
- 6. Add I ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

7. Transfer all of the mixture to a Midi column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of the mixture to a Midi column, be careful not to moisten the rim of a Midi column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

8. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column.

9. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the column membrane has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the Midi column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

- 10. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg (5,000 rpm), follow this;
  - → Discard the filtrate, wipe off any spillage from the thread of the 15 ml tube, and re-insert the Midi column back into the 15 ml tube. Apply 1 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the Midi column and incubate it at room temperature for 15 min.

Insufficient centrifugal force will bring on residual ethanol in Midi column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

II. Place the Midi column into a new 15 ml centrifugation tube (provided). Pipet 300 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge for 5 min at 4,500 xg (5,000 rpm).

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in Midi column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the Midi column membrane for optimal elution of DNA.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

If the weight or the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 100 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problems by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water is higher than 7.0.

12. For higher concentrated yield, re-load the eluate from step 11 into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

For higher total yield, add 300 ul of fresh Buffer AE or distilled water again into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4.500 xg (5.000 rpm)

The first and second eluates can be combined or collected seperately as necessity.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

## PROTOCOL for 100 ~ 250 mg animal tissue

#### Before experiment





Prepare the water bath to 56°C and 70°C

Prepare absolute ethanol

If a precipitate has formed in buffer CL or BL, heat to dissolve completely

1. Homogenize  $100 \sim 250$  mg of tissue as described in step IA, IB or IC, depending on the sample type.

It is most important to weigh the sample accurately.

If the sample is spleen tissue, up to 100 mg can be processed.

Well-disrupted sample will accelerate lysis and decrease the lysis time.

- IA. For soft tissue, such as liver or brain, put up to 250 mg of the tissue into homogenizer, add I ml of distilled water, homogenize thoroughly on ice. Transfer the homogenate into 50 ml conical tube, add 3 ml of Buffer CL, and vortex vigorously to homogenate well.
- IB. If a homogenizer is not available or the tissue is not soft, grind the tissue to a fine powder with liquid nitrogen in a pre-cooled mortar and pestle. Put up to 250 mg of the powdered tissue into 50 ml conical tube. Add 4 ml of Buffer CL and vortex for 30 sec to homogenate completely.
- I C. If neither IA nor IB is available, mince the tissue with sterile sharp blade or scalpel as small as possible. Put up to 250 mg of the tissue into a 50 ml conical tube. Add 4 ml of Buffer CL and pulse-vortex for 30 sec.
  - \*\* Alternatively, tissue sample can be effectively disrupted using some instruments, such as rotor-stator homogenizer or a bead-beater.

    When use these, follow the manufacture s instruction manual.

2. Add 200 ul of Proteinase K solution to the tube. Mix completely by vortexing or pipetting. Incubate at 56°C until the sample is completely lysed.

It is essential to mix the components completely for efficient lysis.

Lysis time varies from 10 min to 3 hr usually depending on the type of tissue and the disruption method. The lysate should become translucent without any particles after complete lysis. Overnight lysis does not influence the preparation. If the sample is lysed in water bath or heating block, vortex occasionally (2  $\sim$  3 times per hour) during incubation to lyse readily. Lysis in shaking water bath, shaking incubator, or agitator would be best for efficient lysis.

3. (Optional:) If RNA-free DNA is required, cool the lysate to room temperature, add 100 ul of RNase solution (100 mg/ml, Cat.No.117-961), vortex to mix thoroughly, and incubate for 3 min at room temperature.

Unless RNase is treated, both DNA and RNA will be co-purified. RNA can inhibit some downstream enzymatic reactions, but not PCR itself.

4. Add 5 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely.

- 5. Incubate at 70°C for 10 min.
- 6. Add 4 ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

7. Transfer all of the mixture to a MAXI column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

Do NOT place the MAXI columns in tilted or bottom-up position even if caps are closed.

8. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply 7 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 50 ml tube before reinserting the MAXI column.

9. Discard the filtrate, re-insert the MAXI column back into the 50 ml tube. Apply 10 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the MAXI column has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the MAXI column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

### 10. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

 $\rightarrow$  Discard the filtrate, wipe off any spillage from the thread of the 50 ml tube, and re-insert the MAXI column back into the 50 ml tube. Apply 3 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the MAXI column and incubate it for 20 min at room temperature.

Insufficient centrifugal force will bring on residual ethanol in MAXI column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

II. Place the MAXI column into a new 50 ml centrifugation tube (provided). Pipet 600 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in MAXI column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the MAXI column membrane for optimal elution of DNA.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yields.

If the weight or the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 300 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA (0.5mM) included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water (at least 7.0) before elution.

12. For higher concentrated yield, re-load the eluate from step 11 into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

For higher total yield, add 600 ul of fresh buffer AE or distilled water again into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

The first and second eluates can be combined or collected seperately as necessity.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

## PROTOCOL for $\sim 1 \times 10^{10}$ gram (-) bacteria

#### Before experiment



Prepare the water bath to 56°C and 70°C

Prepare absolute ethanol

If a precipitate has formed in buffer CL or BL, heat to dissolve completely

I. Harvest bacterial cells (up to 1 x 10<sup>10</sup>) in a 15 ml conical tube by centrifugation at 10,000 xg for 5 min. Discard the supernatant as much as possible.

The number of cells in a bacterial culture varies depending on each strain. When  $A_{600}=I$ , 10 ml of bacterial culture may correspond to  $I\sim 2\times 10^{10}$  cells approximately.

- 2. Resuspend the cell pellet with the residual liquid by flickering or vortexing.
- 3. Add I ml of buffer CL and resuspend completely by pipetting or vortexing.
- 4. Pipet 50 ul of Proteinase K solution (20 mg/ml, provided) and mix completely by pipetting or vortexing.
- 5. Incubate 56°C for 20 min.

Vortex the lysate occasionally during incubation to accelerate lysis. Longer incubation will not affect DNA recovery.

6. (Optional:) If RNA-free DNA is required, cool the mixture to room temperature, add 20 ul of RNase solution (100 mg/ml, Cat.No.117-961), vortex to mix thoroughly, and incubate for 5 min at room temperature.

Unless RNase is treated, both DNA and RNA will be co-purified. RNA can inhibit some downstream enzymatic reactions, but not PCR itself.

7. Add I.2 ml of Buffer BL to the tube. Vortex the tube for 15 sec to mix thoroughly.

For efficient lysis and consistent result, it is essential to mix the sample completely.

- 8. Incubate at 70°C for 10 min.
- 9. Add I ml of absolute ethanol (not provided) to the sample, vortex to mix the sample thoroughly.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

10. Transfer all of the mixture to a Midi column carefully, close the cap, centrifuge for 3 min at 2,000 xg (3,000 rpm).

While transfer of the mixture to a Midi column, be careful not to moisten the rim of a Midi column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

11. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of Buffer BW and centrifuge for 3 min at 2,000 xg (3.000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column.

12. Discard the filtrate and re-insert the Midi column back into the 15 ml tube. Apply 3 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5.000 rpm).

If the column membrane has residual ethanol(originated from buffer TW) associated with it after centrifugation, incubate the Midi column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

### 13. (Optional:) If the centrifugal force applied at previous step is less than $4,500 \times g$ , follow this;

→ Discard the filtrate, wipe off any spillage from the thread of the 15 ml tube, and re-insert the Midi column back into the 15 ml tube. Apply 1 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the Midi column and incubate it at room temperature for 15 min.

Insufficient centrifugal force will bring on residual ethanol in Midi column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

# 14. Place the Midi column into a new 15 ml centrifugation tube (provided). Pipet 300 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge for 5 min at 4,500 xg (5,000 rpm).

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in Midi column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the Midi column membrane for optimal elution of DNA.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

If the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 100 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problems by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water is higher than 7.0.

15. For higher concentrated yield, re-load the eluate from step 14 into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm).

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

For higher total yield, add 300 ul of fresh Buffer AE or distilled water again into the Midi column, close the cap, incubate 5 min at room temperature, and centrifuge for 5 min at 4,500 xg (5,000 rpm)

The first and second eluates can be combined or collected seperately as necessity.

Less than 300 ul of eluate will be obtained from 300 ul of elution buffer, but this has no influence on DNA yields.

## PROTOCOL for $\sim 5 \times 10^{10}$ gram (-) bacteria

#### Before experiment





Prepare the water bath to 56°C and 70°C

Prepare absolute ethanol

If a precipitate has formed in buffer CL or BL, heat to dissolve completely

I. Harvest bacterial cells (up to 5 x 10<sup>10</sup>) in a 50 ml conical tube by centrifugation at 10,000 xg for 5 min. Discard the supernatant as much as possible.

The number of cells in a bacterial culture varies depending on each strain. When  $A_{600}=I$ , I0 ml of bacterial culture may correspond to  $I\sim 2\times I0^{10}$  cells approximately.

- 2. Resuspend the cell pellet with the residual liquid by flickering or vortexing.
- 3. Add 5 ml of buffer CL and resuspend completely by pipetting or vortexing.
- 4. Pipet 200 ul of Proteinase K solution (20 mg/ml, provided) and mix completely by pipetting or vortexing.
- 5. Incubate 56°C for 30 min.

Vortex the lysate occasionally during incubation to accelerate lysis. Longer incubation will not affect DNA recovery.

6. (Optional:) If RNA-free DNA is required, cool the mixture to room temperature, add 50 ul of RNase solution (100 mg/ml, Cat.No.117-961), vortex to mix thoroughly, and incubate for 5 min at room temperature.

Unless RNase is treated, both DNA and RNA will be co-purified. RNA can inhibit some downstream enzymatic reactions, but not PCR itself.

- 7. Add 6 ml of buffer BL to the tube. Vortex the tube to mix thoroughly.
- 8. Incubate at 70°C for 10 min.
- 9. Add 5 ml of absolute ethanol to the sample, mix thoroughly by inverting or pulse-vortexing.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

10. Transfer all of the mixture to a MAXI column carefully, close the cap, centrifuge for 3 min at 2,000xg (3,000 rpm).

While transfer of mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

Do NOT place the MAXI columns in tilted or bottom-up position even if caps are closed.

II. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply 7 ml of Buffer BW and centrifuge for 2 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column.

12. Discard the filtrate, re-insert the MAXI column back into the 50 ml tube. Apply 10 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the MAXI column has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the MAXI column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

### 13. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

→ Discard the filtrate, wipe off any spillage from the thread of the 50 ml tube, and re-insert the MAXI column back into the 50 ml tube. Apply 3 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the MAXI column and incubate it for 20 min at room temperature.

Insufficient centrifugal force will bring on residual ethanol in MAXI column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

# 14. Place the MAXI column into a new 50 ml centrifugation tube (provided). Pipet 600 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in MAXI column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the MAXI column membrane for optimal elution of DNA.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yields.

If the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 300 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA (0.5mM) included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water (at least 7.0) before elution.

15. For higher concentrated yield, re-load the eluate from step 14 into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

For higher total yield, add 600 ul of fresh buffer AE or distilled water again into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

The first and second eluates can be combined or collected seperately as necessity.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

## PROTOCOL for $\sim 5 \times 10^{10}$ gram (+) bacteria

#### Before experiment



Prepare the water bath to 37°C, 56°C and 70°C

Prepare Lysozyme or Lysostaphin

Prepare absolute ethanol

If a precipitate has formed in buffer CL or BL, heat to dissolve completely

All centrifugation should be performed at room temperature

#### **Prepare Enzyme Mixture**

Dissolve the appropriate enzyme (not provided, listed below) with buffer GP just before use. Enzyme mixture should be stored at -20°C (or below) as small aliquots; Ideally, once per an aliquot. Thawed aliquot should be discarded after use.

30 mg/mL lysozyme (LYS702, Bioshop, Canada or equivalents) or/and

300 ug/mL lysostaphin (L7386, Sigma, USA or equivalents)

\* For certain species, such as Staphylococcus, treatment of lysostaphin (final conc. = 300 ug/mL) may be required for efficient lysis instead of (or with) lysozyme. However, lysozyme is sufficient to lyse the cell wall for most gram positive bacterial strains.

1. Harvest bacterial cells (up to 5 x 10<sup>10</sup>) in a 50 ml conical tube by centrifugation at 10,000 xg for 5 min. Discard the supernatant as much as possible.

The number of cells in a bacterial culture varies depending on each strain. When  $A_{600} = 1.10$  ml of bacterial culture may correspond to  $1 \sim 2 \times 10^{10}$  cells approximately.

- 2. Resuspend the cell pellet with the residual liquid by flickering or vortexing.
- 3. Add 5 ml of the prepared enzyme mixture and resuspend completely by pipetting or vortexing.
- 4. Incubate at 37°C for 40 min.

The purpose of this treatment is to weaken the cell wall so that cell lysis can be efficiently taken place.

5. (Optional:) If RNA-free DNA is required, add 50 ul of RNase solution (100 mg/ml, cat.no.117-961), vortex to mix thoroughly, and incubate for 5 min at room temperature.

Unless RNase A is treated, RNA will co-purified with DNA, RNA may inhibit some downstream enzymatic reactions, but not PCR itself.

6. Add 6 ml of buffer BL and 200 ul of Proteinase K solution (20 mg/ml, provided) and vortex vigorously to mix completely.

Incubate 56°C for 30 min and then for a further 30 min at 70°C. Cool to room temperature.

If pathogen, it is strongly recommended to substitute incubating at 70°C for 30 min by incubating at 95°C for 15 min.

### 8. Add 5 ml of absolute ethanol to the sample, mix thoroughly by inverting or pulse-vortexing.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

### 9. Transfer all of the mixture to a MAXI column carefully, close the cap, centrifuge for 3 min at 2,000xg (3,000 rpm).

While transfer of mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

Do NOT place the MAXI columns in tilted or bottom-up position even if caps are closed.

#### I 0. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply 7 ml of buffer BW and centrifuge for 2 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column

## II. Discard the filtrate, re-insert the MAXI column back into the 50 ml tube. Apply 10 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the MAXI column has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the MAXI column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

#### 12. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

→ Discard the filtrate, wipe off any spillage from the thread of the 50 ml tube, and re-insert the MAXI column back into the 50 ml tube. Apply 3 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the MAXI column and incubate it for 20 min at room temperature.

Insufficient centrifugal force will bring on residual ethanol in MAXI column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

#### 13. Place the MAXI column into a new 50 ml centrifugation tube (provided). Pipet 600 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in MAXI column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the MAXI column membrane for optimal elution of DNA.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yields.

If the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 300 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA (0.5mM) included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water (at least 7.0) before elution.

14. For higher concentrated yield, re-load the eluate from step 13 into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

For higher total yield, add 600 ul of fresh buffer AE or distilled water again into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

The first and second eluates can be combined or collected seperately as necessity.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

## PROTOCOL for $\sim 5 \times 10^8$ yeast

#### Before experiment



Prepare the water bath to 37°C, 56°C and 70°C

Prepare Lyticase or Zymolase

Prepare absolute ethanol

1. Add 30  $\sim$  50 ml (up to 5 x 10<sup>8</sup>) of a yeast culture grown in YPD broth to a 50 ml conical tube.

When the value of A600 reaches to 1.0 (generally, log-phase), 30 ml of culture may yield  $100 \sim 300$  ug DNA approximately.

- 2. Harvest cells in a 50 ml conical tube by centrifugation at 5,000 xg for 10 min. Discard the supernatant as much as possible.
- 3. Resuspend the cell pellet with the residual liquid by flickering or vortexing.
- 4. Add 5 ml of buffer YL and resuspend completely by pipetting or vortexing.
- 5. Add 2,000 U of lyticase or 200 U of zymolase and gently pipet to mix completely.

Unit/mg of enzyme will vary depending on the manufactures.

6. Incubate at 37°C for I hr to digest the cell wall.

Incubated cells turn to spheroplasts at this step, and this makes easy the followed lysis step.

7. Centrifuge at 5,000 xg for 10 min and remove the supernatant as much as possible.

- 8. Resuspend the cell pellet with the residual liquid by flickering or vortexing.
- 9. Add 5 ml of buffer CL and resuspend completely by pipetting or vortexing.
- I 0. Pipet 200 ul of Proteinase K solution (20 mg/ml, provided) and mix completely by pipetting or vortexing.
- II. Incubate 56°C for 30 min.

Vortex the lysate occasionally during incubation to accelerate lysis. Longer incubation will not affect DNA recovery.

12. (Optional:) If RNA-free DNA is required, cool the mixture to room temperature, add 50 ul of RNase solution (100 mg/ml, Cat.No.117-961), vortex to mix thoroughly, and incubate for 5 min at room temperature.

Unless RNase is treated, both DNA and RNA will be co-purified. RNA can inhibit some downstream enzymatic reactions, but not PCR itself.

- 13. Add 6 ml of buffer BL to the tube. Vortex the tube to mix thoroughly.
- 14. Incubate at 70°C for 10 min.
- 15. Add 5 ml of absolute ethanol to the sample, mix thoroughly by inverting or pulse-vortexing.

It is essential to mix the sample completely for efficient binding.

A white thread-like strands can be formed in the lysate. It is essential to transfer all of the lysate including this to a Midi column at next step.

#### 16. Transfer all of the mixture to a MAXI column carefully, close the cap, centrifuge for 3 min at 2,000xg (3,000 rpm).

While transfer of mixture to the MAXI column, be careful not to moisten the rim of MAXI column.

If the mixture has not passed completely through the membrane, centrifuge again at higher speed until all of the solution has passed through.

Do NOT place the MAXI columns in tilted or bottom-up position even if caps are closed

#### 17. Discard the filtrate and re-insert the MAXI column back into the 50 ml tube. Apply 7 ml of Buffer BW and centrifuge for 2 min at 2,000 xg (3,000 rpm).

Wipe off any spillage from the thread of the 15 ml tube before reinserting the Midi column

#### 18. Discard the filtrate, re-insert the MAXI column back into the 50 ml tube. Apply 10 ml of buffer TW, and centrifuge for 15 min at 4,500 xg (5,000 rpm).

If the MAXI column has residual ethanol (originated from buffer TW) associated with it after centrifugation, incubate the MAXI column for 15 min at room temperature to evaporate residual ethanol.

This residual ethanol can decrease DNA yield significantly and it also can inhibit some downstream reactions.

## 19. (Optional:) If the centrifugal force applied at previous step is less than 4,500 xg, follow these;

→ Discard the filtrate, wipe off any spillage from the thread of the 50 ml tube, and re-insert the MAXI column back into the 50 ml tube. Apply 3 ml of absolute ethanol and centrifuge for additional 15 min at available full speed. Remove the MAXI column and incubate it for 20 min at room temperature.

Insufficient centrifugal force will bring on residual ethanol in MAXI column membrane, followed by poor DNA recovery. At least, 4,000 xg is required for proper DNA recovery.

# 20. Place the MAXI column into a new 50 ml centrifugation tube (provided). Pipet 600 ul of buffer AE or distilled water onto a center of membrane and close the cap. Incubate for 5 min at room temperature. Centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Before this elution step, it is strongly recommended that any residual ethanol originated from buffer TW should not remain in MAXI column membrane. It can be checked by smell. Residual ethanol disturbs proper DNA preparation as follows; poor DNA recovery, low purity, inhibition of subsequent reactions and etc.

Ensure that the Buffer AE or distilled water is dispensed directly onto the center of the MAXI column membrane for optimal elution of DNA.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yields.

If the cell number of starting sample is small, less volume of buffer AE or distilled water can be applied. However, do not reduce the elution volume below 300 ul.

For long-term storage, eluting in buffer AE is recommended. But, EDTA (0.5mM) included in the buffer AE may inhibit subsequent enzymatic reactions, so you can avoid such latent problem by using distilled water (>pH7.0) or Tris-Cl (>pH8.5). When using water for elution, make sure the pH of water (at least 7.0) before elution.

21. For higher concentrated yield, re-load the eluate from step 20 into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

For higher total yield, add 600 ul of fresh buffer AE or distilled water again into the MAXI column, close the cap, incubate 5 min at room temperature, and centrifuge at 4,500 xg (5,000 rpm) for 5 min.

The first and second eluates can be combined or collected seperately as necessity.

Less than 600 ul of eluate will be obtained from 600 ul of elution buffer, but this has no influence on DNA yield.

#### **Troubleshooting Guide**

Facts	Possible Causes	Suggestions
Low or no recovery	Low cells in the sample	Some blood may have low concentration of white blood cells. Increase the sample volumes and load the column for binding several times. Reduce the elution volume to minimum. If possible, draw new blood sample and repeat the DNA purification with new sample.
	Too much starting sample	If too much cells present in the sample, reduce the starting sample weight or volume, or increase the volume of buffers to double.
	Starting material was too old or mis-stored.	Best result is obtained when using a fresh sample. DNA yield is dependent on the type, size, age and storage condition of starting material. Lower yield may be obtained from the material that has not been appropriately stored. If blood, samples that have been stored at 4°C for more than 5 days may give reduced yield.
	Inefficient lysis	Inefficient lysis may be due to insufficient mixing with Buffer BL, too much cells in the starting sample, or degenerated Proteinase K. After addition of Buffer CL or BL in procedures, vortex the mixture vigorously and immediately to mix completely.
	Weaken activity of Proteinase K caused by mis-storage or out- of-date	Proteinase K should be stored at 4°C for maintenance of proper activity. However, it is recommended to store at -20°C for prolonged preservation of its activity. Lysis can not be done properly with degenerated Proteinase K.

Facts	Possible Causes	Suggestions
Low or no recovery	The centrifugation step in procedure is performed on the fixed-angle centrifuge.	The centrifugation step in procedure must be carried out on swinging-bucket-have centrifuge. If not, DNA yields may be significantly reduced.
	G-force was not reach to 4,500 x g	For proper DNA purification, centrifugal g-force should be reached to $4,500 \times g$ at least. If g-force is under $4,500 \times g$ , refer to the step II.
	Column not incubated for 5 min	After addition of Buffer AE or distilled water, the column should be incubated at room temperature for 5 min before centrifugation.
	Improper eluent	As user's requirement, elution buffer other than Buffer AE can be used. However, the optimal conditions for elution should be considered carefully whether low salt concentration with alkaline pH (7 <ph<9). as="" buffer="" conditions.<="" eluent,="" employed="" ensure="" or="" other="" th="" that="" was="" water="" when=""></ph<9).>
SV column has colored residue associated	Inefficient lysis	Inefficient lysis may cause that colored residue remains on the column membrane. Repeat the procedure after consideration of 'Inefficient lysis' at 'Low or no recovery'.
with it after wash, resulting in colored residue	Incomplete removal of hemoglobin	In case of DNA preparation from certain animal blood, it is hard to remove hemoglobin from their blood. Carry out additional wash step with Buffer BW.

Facts	Possible Causes	Suggestions
Low A <sub>260/280</sub> ratio	Insufficient lysis	Insufficient lysis cause low DNA purity, and is due to incomplete mixing with Buffer BL or ethanol, too much cells in the starting sample, or degenerated Proteinase K. Check these out in next preparations.
	Residual ethanol from Buffer TW remains in eluate	If centrifugation has been performed on fixed-angle rotor or centrifugal force does not reach 4,000 xg, residual ethanol will remain in the column membrane and it causes low purity. To remove residual ethanol completely, refer to the optional step.
	Incomplete removal of hemoglobin	In case of DNA preparation from certain animal blood, it is hard to remove hemoglobin from their blood. Carry out additional wash step with Buffer BW.
High A <sub>260/280</sub> ratio	RNA contamination	RNA may inhibit some downstream enzymatic reactions. If RNA-free DNA is required, carry out the optional step.
Low concentration of DNA in eluate	Low cells in starting sample (too high elution volume)	Increase the volume of starting sample, or reduce the elution volume to minimum or do re-elution with eluate.  Check the last step in procedures for higher concentrations of DNA.
Column clogging	Inefficient lysis	Inefficient lysis may lead to column clogging. About inefficient lysis, check 'Inefficient lysis' at 'Low or no recovery'.
Degraded DNA	Starting sample is too old or mis-stored	Too old or mis-stored sample often yield degraded DNA. Use fresh sample.

Facts	Possible Causes	Suggestions
DNA floats out of well while loading of agarose gel	Residual ethanol from Buffer TW remains in eluate	If centrifugation has been performed on fixed-angle rotor or centrifugal force does not reach 4,500 xg, residual ethanol will remain in the column membrane. Residual ethanol interferes with DNA in sinking into a well because of its low specific gravity. To remove residual ethanol completely, refer to the optional step.
Enzymatic	Low purity of DNA	Check 'Low A <sub>260/280</sub> ratio'
reaction is not performed well with	RNA contamination	RNA may inhibit some downstream enzymatic reactions. If RNA-free DNA is required, carry out the optional step.
purified DNA	Incomplete removal of hemoglobin	In case of DNA preparation from certain animal blood, it is hard to remove hemoglobin from their blood. Carry out additional wash step with Buffer BW.
	High salt concentration in eluate	Ensure that all washing steps were performed just in accordance with the protocols.  Optionally, additional washing with buffer TW can help remove some more salts.
Precipitate in some buffer	Buffer stored in cool ambient condition	For proper DNA purification, any precipitate in a Buffer should be dissolved by incubating it at 37°C (or above) until it disappears.

### **Ordering Information**

Products	Scale	Size	Cat. No.	Туре	Products	Scale	Size	Cat. No.	Туре
	. aTM s					TM .			
	<b>d-Q</b> '''' fo			plasmid DNA	GeneAll® Exgen	for is			
Plasmid Rapidprep mini	50	100-150	spin		mini	100	105-101	spin /	
	200	100-102				250	105-152	vacuum	
GeneAll® <i>Expr</i> e	_TM			DAIA	Blood SV	Midi	26	105-226	spin / vacuum
JeneAll Expre	<b>p</b> jor pi		101-150	DINA			100	105-201	
	50 200	101-150	spin /		MAXI	26	105-316	spin / vacuum	
	mini	1,000	101-102	vacuum			100	106-101	spin /
Plasmid SV	-	26	101-111			mini	250	106-152	vacuum
	Midi	50	101-250	spin /	Cell SV		10	106-310	spin /
		100	101-201	vacuum		MAXI	26	106-326	vacuum
							100	108-101	spin /
GeneAll® <i>Exfec</i> t	tion™					mini	250	108-152	vacuum
for prep	aration of	transfect	ion-grade pla		N 41 11	26	108-226	spin /	
	mini	50	111-150	spin /	Midi	100	108-201	vacuum	
Plasmid LE		200	111-102	vacuum	/	MAXI	10	108-310	spin /
(Low Endotoxin)	Midi	26	111-226	spin /		MAXI	26	108-326	vacuum
		100	111-201	vacuum	vacuum Genomic DNA micro	Ö	50	118-050	spin
Plasmid EF Midi	20	121-220	- snin			100	117-101	spin /	
/E / E /	i iidi								
(Endotoxin Free)	1 1101	100	121-201	Spii i		mini	250	117-152	vacuum
,					Plant SV		250 26	117-152 117-226	vacuum spin /
(Endotoxin Free) GeneAll® <i>Expin</i>		ification	of fragment D	DNA	Plant SV	Midi			
,		ification	of fragment D 102-150	DNA spin /	Plant SV	Midi	26	117-226   117-201   117-310	spin /
GeneAll® <i>Expin</i>	<b>тм</b> for pur	50 200	of fragment D 102-150 102-102	oNA spin / vacuum			26 100 10 26	117-226   117-201   117-310   117-326	spin / vacuum
GeneAll® <i>Expin</i>	<b>тм</b> for pur	50 200 50	of fragment C 102-150 102-102 103-150	spin / vacuum spin /	Soil DNA mini	Midi	26 100 10 26 50	117-226 117-201 117-310 117-326 114-150	spin / vacuum spin /
GeneAll® Expin Gel SV	for pur mini	50 200 50 200	102-150 102-102 103-150 103-102	spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini	Midi  MAXI  mini  mini	26 100 10 26 50	117-226 117-201 117-310 117-326 114-150 115-150	spin / vacuum spin / vacuum spin spin
GeneAll® Expin Gel SV	for pur mini	50 200 50 200 50 200	102-150 102-102 103-150 103-102 113-150	spin / vacuum spin / vacuum spin /	Soil DNA mini	Midi  MAXI  mini	26 100 10 26 50	117-226 117-201 117-310 117-326 114-150	spin / vacuum spin / vacuum spin
GeneAll® Expin  Gel SV  PCR SV	for pur mini mini	50 200 50 200 50 200 50 200	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102	spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA	Midi  MAXI  mini  mini  mini	26 100 10 26 50 50	117-226 117-201 117-310 117-326 114-150 115-150 128-150	spin / vacuum spin / vacuum spin spin spin
GeneAll® Expin  Gel SV  PCR SV	for pur mini mini	50 200 50 200 50 200 50 200	102-150 102-102 103-150 103-150 103-102 113-150 113-102 112-150	spin / vacuum	Soil DNA mini Stool DNA mini	Midi  MAXI  mini  mini  mini	26 100 26 50 50 50	117-226 117-201 117-310 117-326 114-150 115-150 128-150	spin / vacuum spin / vacuum spin spin spin
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV	for pur mini mini mini	50 200 50 200 50 200 50 200	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102	spin / vacuum spin / vacuum spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA GeneAll® GenEx	Midi  MAXI  mini  mini  mini	26 100 10 26 50 50 50	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit	spin / vacuum spin / vacuum spin spin spin
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP	for pur mini mini mini mini	50 200 50 200 50 200 50 200 50 200	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA	Midi  MAXI  mini  mini  mini  mini  TM for isol	26 100 10 26 50 50 50 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105	spin / vacuum spin / vacuum spin spin spin spin spin spin spin
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV	for pur mini mini mini mini	50 200 50 200 50 200 50 200 50 200 0lation o	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA GeneAll® GenEx	MAXI  mini mini mini mini	26 100 10 26 50 50 50 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301	spin / vacuum spin / vacuum spin spin spin spin spin spin spin
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP	for pur mini mini mini mini	50 200 50 200 50 200 50 200 50 200 olation o	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx <sup>TM</sup> Blood	Midi  MAXI  mini  mini  mini  mini  TM for isol	26 100 10 26 50 50 50 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101	spin / vacuum spin / vacuum spin spin spin spin spin spin spin
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP	mini mini mini mini mini	50 200 50 200 50 200 50 200 50 200 0lation of 100 250	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA GeneAll® GenEx	Midi  MAXI  mini  mini  mini  mini  for isol  Sx  Lx  Sx	26 100 26 50 50 50 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105	spin / vacuum spin / vacuum spin / spin spin spin spin spin spin spin solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP	mini mini mini mini mini	50 200 50 200 50 200 50 200 50 200 100 250 26	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx <sup>TM</sup> Blood	Midi  MAXI  mini mini mini  TM for isol  Sx  Lx	26 100 10 26 50 50 6dtion of 100 500 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301	spin / vacuum spin / vacuum spin / spin spin spin spin spin spin spin solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP  GeneAll® Exgen	for pur mini mini mini mini mini mini mini min	50 200 50 200 50 200 50 200 50 200 0lation of 100 250 26 100	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx™ Blood  GenEx™ Cell	Midi  MAXI  mini  mini  mini  mini  for isol  Sx  Lx  Sx	26 100 10 26 50 50 6dtion of 100 500 100 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101	spin / vacuum spin / vacuum spin / spin spin spin spin spin spin spin solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP  GeneAll® Exgen	for pur mini mini mini mini mini mini mini min	50 200 50 200 50 200 50 200 50 200 0lation of 100 250 26 100	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx <sup>TM</sup> Blood	Midi  MAXI  mini mini mini  TM for isol  Sx  Lx  Sx  Lx  Sx	26 100 10 26 50 50 50 action of 100 500 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin / vacuum spin / vacuum spin / spin spin spin spin spin solution solution solution solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP  GeneAll® Exgen	for pur mini mini mini mini mini mini mini Midi	50 200 50 200 50 200 50 200 50 200 50 200 250 25	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx™ Blood  GenEx™ Cell	Midi  MAXI  mini mini mini  for isol  Sx  Lx  Sx  Lx	26 100 10 26 50 50 6dtion of 100 500 100 100 100	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101	spin / vacuum spin / vacuum spin / spin spin spin spin spin solution solution solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP  GeneAll® Exgen	for pur mini mini mini mini mini mini mini Midi	50 200 50 200 50 200 50 200 50 200 60 61 60 250 26 100 10 26 100	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx™ Blood  GenEx™ Cell	Midi  MAXI  mini mini mini  TM for isol  Sx  Lx  Sx  Lx  Sx	26 100 10 26 50 50 50 action of 100 500 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin / vacuum spin / vacuum spin spin spin spin spin spin solution solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP  GeneAll® Exgen	mini mini mini mini mini mini mini mini	50 200 50 200 50 200 50 200 50 200 0lation of 250 26 100 26 100 250	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx™ Blood  GenEx™ Cell	Midi  MAXI  mini mini mini  TM for isol  Sx  Lx  Sx  Lx  Sx	26 100 10 26 50 50 50 action of 100 500 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin / vacuum spin / vacuum spin / spin spin spin spin spin solution solution solution solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP  GeneAll® Exgen	mini mini mini mini mini mini mini mini	50 200 50 200 50 200 50 200 50 200 0lation of 100 250 26 100 10 26 100 250 26	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx™ Blood  GenEx™ Cell	Midi  MAXI  mini mini mini  TM for isol  Sx  Lx  Sx  Lx  Sx	26 100 10 26 50 50 50 action of 100 500 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin / vacuum spin / vacuum spin / spin spin spin spin spin solution solution solution solution solution
GeneAll® Expin  Gel SV  PCR SV  CleanUp SV  Combo GP  GeneAll® Exgen  Tissue SV	mini mini mini mini mini mini mini mini	50 200 50 200 50 200 50 200 50 200 0lation of 250 26 100 26 100 250	of fragment D 102-150 102-102 103-150 103-102 113-150 113-102 112-150 112-102 f total DNA 104-101 104-152 104-226 104-201 104-310 104-326 109-101 109-152	spin / vacuum	Soil DNA mini Stool DNA mini Viral DNA / RNA  GeneAll® GenEx  GenEx™ Blood  GenEx™ Cell	Midi  MAXI  mini mini mini  TM for isol  Sx  Lx  Sx  Lx  Sx	26 100 10 26 50 50 50 action of 100 500 100 100 100 500	117-226 117-201 117-310 117-326 114-150 115-150 128-150 total DNA wit 220-101 220-105 220-301 221-101 221-105 221-301 222-101 222-105	spin / vacuum spin / vacuum spin / spin spin spin spin spin solution solution solution solution solution

Products	Scale	Size	Cat. No.	Туре	Products	Scale	Size	Cat. No.	Туре
GeneAll® <b>GenE</b> x	τм <sub>for is</sub>	olation of	total DNA		GeneAll® <b>Amp</b> O	NE <sup>TM</sup> fo	r PCR aı	mblification	
Bellevill Consul	Sx	100	227-101		ochera - m-p -		250 L		
GenEx <sup>™</sup> Plant	Mx	100	227-201	solution	Tag DNA polymeras	Tag DNA polymerase			(2.5 U/µℓ)
-	Lx	100	227-301		, , ,		1,000 L		-(=:= =, µ==)
	Sx	100	228-101			-		J 502-025	
GenEx <sup>™</sup> Plant plus!	Mx	50	228-250	solution	α-Tag DNA polyme	lpha-Tag DNA polymerase			_ (2.5 U/ <b>µℓ</b> )
	Lx	20	228-320				1,000 L		(/
GeneAll® <i>DirEx</i> ™	<sup>M</sup> serie	s					250 L		
. , ,	ration of	PCR-ten	nplate without	extraction	α-Pfu DNA polyme	erase	500 L		- (2.5 U/ <b>ມໃ</b> )
$DirEx^{TM}$		100	250-101	solution	22 / 12 2 / 12 / 12 / 13		1,000 L		(=:= = / [==]
DirEx <sup>™</sup> Fast-Tissue		96 T	260-011	solution			250 L		
DirEx <sup>™</sup> Fast-Culture	d cell	96 T	260-021	solution	Fast-Pfu DNA		500 L		(2.5 U/µℓ)
$DirEx^{TM}$ Fast-Whole by	olood	96 T	260-031	solution	polymerase		1,000 L		(=:= = / [, = /
$DirEx^{TM}$ Fast-Blood st	ain	96 T	260-041	solution			250 L		
DirEx <sup>™</sup> Fast-Hair		96 T	260-051	solution	Hotstart Taq DNA		500 L		(2.5 U/ <b>μℓ</b> )
$DirEx^{TM}\mathit{Fast} ext{-}Buccal\;s$	wab	96 T	260-061	solution	polymerase		1,000 L		
$DirEx^{TM}\mathit{Fast}\text{-}Cigarett$	е	96 T	260-071	solution			20 µl		) lyophilized
							50 µl		
GeneAll® <b>RNA</b> s	eries †	or preper	ation of total	RNA	Taq Premix	96 tube	s 20 µl		
D: E IM		100	301-001	1.2			50 ul		solution
RiboEx <sup>™</sup>	mini	200	301-002	solution			20 µl		
Hybrid-R <sup>™</sup>	mini	100	305-101	spin					- lyophilized
Hybrid-R <sup>™</sup> Blood RN	IA mini	50	315-150	spin	lpha-Taq Premix	96 tube	50 μ <b>l</b> 522-500 ′ ′ 20 μ <b>l</b> 527-200		
Hybrid-R <sup>TM</sup> miRNA	mini	50	325-150	spin			50 μl		solution
		100	302-001				20 µl		
RiboEx <sup>™</sup> LS	mini	200	302-002	solution	HS-Tag Premix	96 tube			solution
Riboclear™	mini	50	303-150	spin	ms-lag rreillix	70 lube			lyophilized
Riboclear <sup>TM</sup> plus!	mini	50	313-150	spin	 α-Pfu Premix	96 tube	20 µl		solution
Ribospin <sup>™</sup>	mini	50	304-150	spin	Tag Premix (w/o dye)	96 tube			lyophilized
Ribospin ™vRD	mini	50	302-150	spin		96 tube			, ,
Ribospin <sup>TM</sup> vRD <i>plus</i>	! mini	50	312-150	spin	dNTPs mix		500 µl	509-020	2.5 mM eac
Ribospin ™vRD II	mini	50	322-150	spin	dNTPs set (set of dATP, dCTP, dGTP ar	nd dTTP)	1 ml x 4 tubes	509-040	100 mM
Ribospin ™ Plant	mini	50	307-150	spin				* Fach dNITI	Ps is available
Ribospin <sup>™</sup> Seed / Fruit	mini	50	317-150	spin				Lacii (IIN I I	s is available

Allspin™

RiboSaver™

mini 50

mini 100

306-150

351-001 solution

spin

Products	Scale	Size	Cat. No.	Туре	
----------	-------	------	----------	------	--

#### GeneAll® AmpMaster<sup>TM</sup> for PCR amplification

Tag Master mix	0.5 ml x 2 tubes	541-010	solution
iaq i laster mix	0.5~mlx10tubes	541-050	solution
α-Taq Master mix	0.5 ml x 2 tubes	542-010	solution
	0.5 ml x 10 tubes	542-050	solution
	0.5 ml x 2 tubes	545-010	solution
HS-Taq Master mix	0.5 ml x 10 tubes	545-050	solution
α-Pfu Master mix	0.5 ml x 2 tubes	543-010	solution
	0.5 ml x 10 tubes	543-050	solution

#### GeneAll® HyperScript<sup>TM</sup> for Reverse Transcription

ase 10,000 U	601-100	solution
$0.5~\mathrm{ml} \times 2~\mathrm{tubes}$	601-710	solution
0.5 ml × 2 tubes	601-730	solution
$0.5~\mathrm{ml} \times 2~\mathrm{tubes}$	601-740	solution
96 tubes, 20 $\mu l$	601-602	solution
96 tubes, 20 μl	601-632	solution
96 tubes, 20 μl	601-642	solution
$0.5~\mathrm{ml} \times 2~\mathrm{tubes}$	602-110	solution
96 tubes, 20 μl	602-102	solution
50 reaction	605-005	solution
10,000 U	605-010	solution
4,000 U	605-004	solution
	0.5 ml × 2 tubes 0.5 ml × 2 tubes 0.5 ml × 2 tubes 96 tubes, 20 μl 96 tubes, 20 μl 0.5 ml × 2 tubes 96 tubes, 20 μl 50 reaction	0.5 ml × 2 tubes 601-710  0.5 ml × 2 tubes 601-730  0.5 ml × 2 tubes 601-740  96 tubes, 20 μl 601-602  96 tubes, 20 μl 601-642  0.5 ml × 2 tubes 602-110  96 tubes, 20 μl 602-102  50 reaction 605-005

#### GeneAll® Protein series

ProtinEx <sup>™</sup> Animal cell / tissue	100 ml	701-001	solution
PAGESTA <sup>™</sup> Reducing 5X SDS-PAGE Sample Buffer	I ml × I0 tubes	751-001	solution

#### Size Cat. No. Products

#### GeneAll® STEADi™ for automatic nucleic acid puritication

STEADi <sup>™</sup> 12 Instrument		GST012
STEADi <sup>™</sup> 24 Instrument		GST024
STEADi <sup>™</sup> Genomic DNA Cell / Tissue Kit	96	401-104
STEADi™ Genomic DNA Blood Kit	96	402-105
STEADi <sup>™</sup> Bacteria DNA Kit	96	403-106
STEADi <sup>™</sup> Total RNA Kit	96	404-304
STEADi <sup>™</sup> Viral DNA/RNA Kit	96	405-322
STEADi <sup>™</sup> CFC Seed DNA/RNA Kit	96	406-C02

#### Note



#### www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro Songpa-gu, Seoul, South Korea 138-859 E-mail: sales@geneall.com

Tel: 82-2-407-0096 Fax: 82-2-407-0779

©2016 GeneAll Biotechnology, All right reserved

Edited by **Bnp**Designed by **KDY**2010.05