

NEW

LED Digital Orbital Shaker

DUAB

GeneAll



Specifications

SK-0180-S

Dimension (W × D × H)	294 × 348 × 123 mm	
Weight	3.1 kg	
Work plate dimension (W × D)	268 × 268 mm	
Max. load capacity (with platform)	2 kg	
Orbital diameter	20 mm	
Display	LED	
Speed range	40~200 rpm	
Time settings range	1 min~99 h 59 min	
Protection class	IP21	
Motor type	Brushless DC motor	
Voltage	100-240 V, 50/60 Hz	
Power	20 W	
Motor input/output power	20 W / 16 W	
Ambient temperature and humidity	5~40°C, 80% RH	
Components	Main unit	1
	Power cable	1
	User Manual	1
	Warranty card	1
	Non-sleep mat	1
	Fixed string	4

✓ 우수한 가성비

₩ 300,000

*부가세 별도

✓ 합리적인 가격



YouTube

* 제품 사용 설명 영상을 확인하세요!

Ordering Information

Cat. No.	Model	Product	Price
8072210201	SK-0180-S	LED Digital Orbital Shaker	300,000원

01 Intuitive LED Display

✓ 시간

미 설정 (00:00) 시 무한 작동
설정 시 최소 1분부터 최대 99시간 59분까지 가능

✓ 속도

최소 40~200 rpm까지 조절 가능

02 Tabletop Size

✓ 294 x 348 mm로 좁은 공간에서도 사용 가능

✓ 100 ø cell culture dish 기준 최대 5개 적재 가능

✓ 다양한 용기 적재 가능 (dish, flask, PCR plate 등)

✓ 최대 적재 무게 : 2 kg


03 Convenient Functions

✓ Brushless DC motor의 조용하고 강력한 shaking

✓ 작동 중 시간 및 속도 조절 가능

✓ IP21 보호등급의 방수 처리

✓ 미끄럼방지 매트 및 고정 스트링 기본 제공



매트



스트링

Memo